# **Growth Inhibition Assay - Summary**

### **Eye Irritation**

The growth inhibition assay determines the inhibition of cell proliferation of established cells lines after exposure to an ocular irritant.

## **Objective & Application**

The growth inhibition assay is a cytotoxicity test to investigate the effect induced by ocular irritants on cell proliferation. The assay has been proposed as an alternative to the Draize rabbit eye irritation test (OECD TG 405, 2012 and Method B.5 of Annex to Commission Regulation 440/2008/EC, EU 2008) to investigate ocular irritancy induced by chemicals (including surfactants, metal salts, alcohols, ketones, and acetates), cosmetic ingredients and consumer products (Reinhardt et al., 1985; Bracher et al., 1987; Reinhardt et al., 1987; Kennah et al., 1989; Horwath-Winter et al., 2004).

## **Basis of the Method**

The effect of an ocular irritant in vitro is assessed by counting the number of cells after treatment. Reductions in cell numbers in comparison to controls are either the result of a reduced cell division rate or/and cell death. These effects are considered indicative for cytotoxic effects on corneal epithelium of the eye.

The growth inhibition assay is used with long-term (48 h) (Reinhard et al., 1985) or with short-term (30 min) exposure periods (Kennah et al., 1989). The inhibition of cell growth is determined directly after the 48 h treatment by Reinhardt et al. (1985), while Kennah et al. (1989) detected the growth inhibition 24 h after treatment. Reinhardt et al. (1985) used fibroblastic cell lines of hamster (BHK) and human origin (Keller cells, MRC-5). Kennah et al. (1989) applied the assay using a murine fibroblast cell line (3T3). The concentration of a test substance resulting in a cell number reduced to 50% of that of controls is called GI-50 (Reinhardt et al., 1987; Kennah et al., 1989).

# **Experimental Description**

### **Biological and Endpoint Measurement:**

CELL NUMBER: is determined with an coulter counter

### **Endpoint Value:**

GI-50: the concentration of a test substance (in mM, org/ml) reducing the cell number to 50% of that of controls

### **Experimental System:**

3T3 FIBROBLASTS (MOUSE): commercially available murine fibroblastic cell line

BHK-21/C13: commercially available fibroblastic cell line derived from Baby Hamster Kidney (BHK)

KELLER CELLS: immortalized human diploid fibroblastic cell line derived from an arm biopsy

MRC-5 CELLS: commercially available human fibroblastic cell line derived from embryonic lung

Cells are seeded in culture plates and incubated for 1 - 2 days prior to exposure. Test substance dilutions are applied to the monolayer for 48 h for long-term exposure (Reinhardt et al., 1985), respectively 30 min for short-term exposure (Kennah et al., 1989). Untreated cultures serve as negative controls. Cells are counted with a coulter counter after trypsination of cells. Reinhardt et al. (1985) determined the lowest concentration of a test substance that induced a significant effect in the cell growth inhibition assay (GIlow, cell numbers at the end of the incubation period in controls and exposed samples were expressed as percentages of the cell number at the beginning of the incubation period, no further information on how these data were transformed to the GIlow was given by the authors).

Kennah et al. (1989) exposed the cells for 30 min and incubated the cells for another 24 h in fresh medium. The authors expressed differences between numbers of treated cells versus those of untreated controls according to the following formula (Kennah et al., 1989):

% GI = (no. of cells of control - no of cells of sample) / (no of cells of control) x 100 The endpoint value GI-50 is determined by linear regression analysis from individual values for % GI, obtained at various test concentrations.

## **Data Analysis/Prediction Model**

Kennah et al. (1989) developed a prediction model for the short-term exposure (30 min) of the growth inhibition assay. The authors performed a linear regression analysis correlating the log GI-50 values with the percentages of corneal swelling *in vivo*. The following equation was derived by Kennah *et al.* (1989):

% corneal swelling = -29.98 (log GI-50) + 135.57

Increasing percentages of corneal swelling in vivo were correlated with increasing severeness of eye irritation by the same authors in a previous paper (Kennah et al., 1989) in the following way:



So far, no prediction model is available for the long-term exposure of the growth inhibition assay.

## **Test Compounds**

Chemicals (including surfactants, metal salts, alcohols, ketones, and acetates), cosmetic ingredients, and consumer products.

## **Modifications**

Horwath-Winter et al. (2004) used primary human conjunctival fibroblasts and exposed these cells for up to 72 h.

### **Discussion**

Using the long-term exposure assay, Reinhardt et al. (1985) evaluated three cytotoxicity assays (cell detachment, cloning efficiency test, and growth inhibition test) using three different cell lines. According to the authors, the ranking order of the cytotoxicity of most chemicals was comparable for all three tests using the different cell lines. The authors reported that the *in vitro* results of the growth inhibition assay correlated well (correlation coefficient was 0.92) with human data (Reinhardt et al., 1985). In another investigation of the same working group three cytotoxicity tests (cell detachment assay, growth inhibition test, membrane permeability assay) and the hen's egg test (HET-CAM, see method summary for more details) were used to classify 24 different tensides (Reinhardt et al., 1987). According to the authors, the growth inhibition assay did not correlate with the *in vivo* data obtained from studies with guinea pigs, which were performed in analogy to the Draize rabbit eye test using guinea pigs as test animals.

Bracher et al. (1987) compared the performance of four in vitro cell toxicity assays (cell detachment assay, the growth inhibition test (long-term exposure), membrane permeability assay, NRU assay), when used with 26 cosmetic ingredients. From the comparison with *in vivo* data Bracher et al. (1987) concluded that the growth inhibition assay produced a large number of false predictions (especially false negatives) and was unable to distinguish minimally irritating substances from strong irritating, respectively from non-irritating ones.

According to Kennah et al. (1989) the short-term exposure growth inhibition assay showed a good correlation between *in vitro* and *in vivo* results for surfactants and alcohols. According to the authors, the correlation between *in vitro* and *in vivo* results was poor for substances of other chemical classes such as ketones, acetates, and aromatics (Kennah et al., 1989).

Selling and Ekwall (1985) developed a MIT to determine the proliferation rate of HeLa cells in the presence of ocular irritants. The authors reported that they could determine cell proliferation, respectively growth inhibition by observing the colour change of the culture medium. After 7 days of culture an increase of cell number (proliferation) is characterized by colour change (indicator phenol red) from red into yellow/orange (pH 6.0 – 6-5), since acid metabolites accumulate in the supernatant of proliferating cells, thus altering the pH of the culture medium. Cultures which do not change the colour of the medium were regarded as inhibited, while a colour change into violet (pH 8.0) was regarded as

complete inhibition of cell proliferation by Selling and Ekwall (1985).

#### **Status**

So far, in the scientific literature, neither evaluation nor validation studies have been reported for the growth inhibition assay.

### **Abbreviations & Definitions**

BHK: Baby Hamster Kidney EC: European Commission EU: European Union GI: Growth Inhibition HET-CAM: Hen's Egg Test on the Chorioallantoic Membrane NRU: Neutral Red Uptake MIT: Metabolic Inhibition Test OECD: Organisation for Economic Co-operation and Development TG: Test Guideline

## **Bibliography**

- Bracher M., Faller C., Spengler J and Reinhardt C.A. (1987) Comparison of in vitro cell toxicity with in vivo eye irritation. Molecular Toxicology 1, 561-570
- Christian, M. S., and Diener, R. M. (1996) Soaps and Detergents - alternatives to animal eye irritation tests. Journal of the American College of Toxicology 15, 1-44
- Horwath-Winter J., Simon M., Kolli H., Trummer G. and Schmut O. (2004) Cytotoxicity evaluation of soft contact lens care solutions on human conjunctival fibroblasts. Ophthalmologica 218(6), 385-389
- Kennah, H. E., Albulescu, D., Hignet, S., and Barrow, C. S. (1989) A critical evaluation of predicting ocular irritancy potential from an in vitro cytotoxicity assay. Fundamental and Applied Toxicology 12, 281-290
- OECD (2012) OECD Test Guideline No. 405: Acute Eye Irritation/Corrosion. First Adopted in 1981. Last update adopted on 2nd October 2012. **[Link to document](http://www.oecd-ilibrary.org/environment/test-no-405-acute-eye-irritation-corrosion_9789264185333-en)** (last access 29.10.2012) **[OECD Guidelines for the Testing of Chemicals, Section 4, Health Effects](http://www.oecd-ilibrary.org/environment/oecd-guidelines-for-the-testing-of-chemicals-section-4-health-effects_20745788)**
- Reinhardt Ch.A., Pelli D.A. and Zbinden G. (1985) Interpretation of cell toxicity data for the estimation of potential irritation. Food and Chemical Toxicology 23, 247-252
- Reinhardt, C. A., Aesbacher, M., Bracher, M., and Spengler, S. (1987) Validation of three cell toxicity tests and the hen's egg test with guinea pig eye and human skin irritation data. Alternative Methods in Toxicology 5, 463-469
- Selling, J. and Ekwall, B. (1985) Screening for eye irritancy using cultured HeLa cells. Xenobiotica 15, 713-717

Last update: February 2010