




European
Commission

STANDARD OPERATING PROCEDURE

T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin, version 1.0

EURL ECVAM validation study of a battery of mechanistic methods relevant for the detection of chemicals that can disrupt the thyroid hormone system

Roszak J.



2023

This Standard Operating Procedure (SOP) has been prepared within the context of a collaboration agreement with the Joint Research Centre (JRC) Directorate for Health, Consumers and Reference Materials (Chemicals Safety and Alternative Methods Unit F3 / EURL ECVAM), for the validation of mechanistic methods to identify potential modulators of thyroid hormone signalling. It aims to provide evidence-based scientific support to the European policymaking process. The contents of this publication do not necessarily reflect the position or opinion of the European Commission. Neither the European Commission nor any person acting on behalf of the Commission is responsible for the use that might be made of this publication. For information on the methodology and quality underlying the data used in this publication for which the source is neither Eurostat nor other Commission services, users should contact the referenced source. The designations employed and the presentation of material on the maps do not imply the expression of any opinion whatsoever on the part of the European Union concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries.

This SOP is part of a series of 3 SOPs used to perform the "T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin":

1. **SOP "T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin" v1.0** (used in Part 1 of the validation study)
2. SOP "Handling, Maintenance and Quality Control of the GH3 cell line" v1.0 (used in Part 1 of the validation study)
3. SOP "Determination of cell proliferation in T screen assay" v1.0 (used in Part 1 of the validation study)

The method was developed by Arno Gutleb, LIST (Luxembourg) and subsequently implemented by the EU-NETVAL test facility NIOM (Poland) within the validation study.

Contact information

Method developer:

Arno Gutleb
PhD, ERT
Luxembourg Institute of Science and Technology (LIST)
Environmental Research and Innovation (ERIN) Department
41, rue du Brill
L-4422 Belvaux
Grand-duchy of Luxembourg
Tel: +352 275 888 481
e-mail: arno.gutleb@list.lu

EU-NETVAL facility

Joanna Roszak
National Center for Alternative Methods to Toxicity Assessment (CMA)
Department of Translational Research
Nofer Institute of Occupational Medicine
91-348 Lodz, 8 Sw. Teresy St, POLAND
e-mail: Joanna.Rozsak@imp.lodz.pl
<https://www.imp.lodz.pl/krajowe-centrum-metod-alternatywnych-do-oceny-toksycznosci>

EU Science Hub

<https://joint-research-centre.ec.europa.eu>

JRC133181

Ispra: European Commission, 2023
© European Union, 2023



The reuse policy of the European Commission documents is implemented by the Commission Decision 2011/833/EU of 12 December 2011 on the reuse of Commission documents (OJ L 330, 14.12.2011, p. 39). Unless otherwise noted, the reuse of this document is authorised under the Creative Commons Attribution 4.0 International (CC BY 4.0) licence (<https://creativecommons.org/licenses/by/4.0/>). This means that reuse is allowed provided appropriate credit is given and any changes are indicated.

For any use or reproduction of photos or other material that is not owned by the European Union, permission must be sought directly from the copyright holders. The European Union does not own the copyright in relation to the following elements:

- Cover page illustration, © BioRender.com

How to cite this report: Roszak J., Standard Operating Procedure for T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin, version 1.0, applied in Part 1 of the EURL ECVAM thyroid validation study, European Commission, Ispra, 2023, JRC133181.

1
2
3
4
5
6
7
8
9
10
11
12
13
14

Standard Operating Procedure

T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin

Author: EU-NETVAL laboratory NIOM, Joanna Roszak

Version: 01

Date: 03.03.2021



CONFIDENTIAL

15 **Contents**

16 1. Description of the method 3

17 1.1 Scope and Limitations of the Application 3

18 1.2 Safety 3

19 1.3 Required Standard Operating Procedures 3

20 1.4 Test System 3

21 1.5 Apparatus and Materials 3

22 1.6 Reagents and Chemicals 4

23 1.7 Reference items, positive and negative control Item(s) 4

24 1.8 Preparations of media and reagents 5

25 1.8.1 Cell culture medium (cDMEM/F12; completed cDMEM/F12) 5

26 1.8.2 PCM medium 5

27 1.8.3 PCM medium +4x EC₅₀ concentration of T3 (PCM-T3 medium) 6

28 1.8.4 Stock and working solutions of chemicals 6

29 1.8.5 Stock solutions of reagents 7

30 2. Verification of the EC₅₀ value of T3 7

31 2.1 Seeding the GH3 cells onto 96-well plate 7

32 2.2 Exposure of the cells to the range of T3 concentrations 7

33 2.3 Measurement of cell proliferation 8

34 2.4 Acceptance criteria 8

35 3. T-Screen experimental procedure 9

36 3.1 Preparation of Test Item 9

37 3.1.1 Stock solutions 9

38 3.1.2 Work solutions to determine solubility 9

39 3.2 Pre-screen experiment 10

40 3.2.1 Seeding the GH3 cells onto 96-well plate 10

41 3.2.2 Preparation of test, reference and control item solutions 10

42 3.2.3 Exposure of the cells 11

43 3.2.4 Measurement of cell proliferation and data analysis 12

44 3.2.5 Acceptance criteria for pre-screen assay 12

45 3.2.6 Identification of the range concentrations of TI for the next experiments 12

46 3.3 T-Screen assay 15

47 3.3.1 Seeding the GH3 cells onto 96-well plate 15

48 3.3.2 Preparation of test, reference and control item solutions 15

49 3.3.3 Exposure of the cells 15

50 3.3.4 Measurement of cell proliferation 20

51 3.4 Data Analysis and Calculations 20

52 3.4.1 Data analysis for AGONISM plate 20

53 3.4.2 Data analysis for ANTAGONISM plate 21

54 3.4.3 Relative potency of test items 23

55 3.5 Acceptance criteria 24

56 3.5.1 Acceptance criteria for AGONISM plate 24

57 3.5.2 Acceptance criteria for ANTAGONISM plate 24

58 4. Reference(s) 24

59

60 *Figure 1.* Scheme of the plate layout for verification of the EC₅₀ value of T3 8

61 *Figure 2.* The T-screen assay workflow with absorbance of AlamarBlue dye detection. 9

62 *Figure 3.* A schema of test item preparation for pre-screen experiment and the T-Screen test. 13

63 *Figure 4.* Plate layout for pre-screen experiment 14

64 *Figure 5.* Plate layout for the 8a T-Screen assay 17

65 *Figure 6.* Exemplary results obtained for (A) AGONISM plate 1 and (B) ANTAGONISM plate 1 23

66 **1. Description of the method**

67 The T-Screen represents an *in vitro* bioassay based on thyroid hormone (TH) dependent cell
68 proliferation of a rat pituitary tumour cell line (GH3) in serum-free medium. It can be used to study
69 interference of compounds with TH at the cellular level, thus bridging the gap between limitations of
70 assays using either isolated molecules (enzymes, transport proteins) or complex *in vivo* experiments
71 with all the complex feedback mechanisms present. Test items are tested both in the absence and
72 presence of TH (EC₅₀ concentration of T3) to test for both agonistic and antagonistic potency.
73 GH3 cell growth is increased in the presence of TH agonists and decreased in the presence of TH plus
74 TH antagonists. Cell growth is measured with AlamarBlue/Resazurin cell proliferation assay using a
75 standard plate reader. In this method a colorimetric assay is used, where resazurin is reduced from a
76 blue oxidized form into its violet reduced form of resorufin. The change of colour can be detected as
77 a change in absorbance using a microplate reader.

78 **1.1 Scope and Limitations of the Application**

79 The T-screen assay is used for testing of compounds for TH receptor (TR)-mediated effects.

80 **1.2 Safety**

81 The use of endocrine disrupting chemicals can be extremely hazardous, and precautions such as
82 using gloves, protective goggles and masks under a laminar flow hood should always be taken while
83 performing chemical treatments.

84 **1.3 Required Standard Operating Procedures**

- 85 • *Handling and Maintenance of GH3 cell line*
- 86 • *Determination of cell proliferation in T-screen assay*

87 **1.4 Test System**

- 88 • GH3-cell line – a rat pituitary tumour cell line, (ATCC® CCL-82.1™)

89

90 The GH3 cells may be used when the following requirements are met:

91 **1.4-1** Test system must be free of microbial and mycoplasma contamination

92 **1.4-2** The doubling time of the GH3 cells in cDMEM/F12 medium should be 42 ± 5h

93 **1.4-3** Cell number in PCM after 72 ± 1 h and 96 ± 1 h of culture should be at least 40% lower
94 that cell number in cDMEM/F12, when determined simultaneously.

95 **1.4-4** Cells should be used in passage from 3 to 20

96 **1.5 Apparatus and Materials**

- 97 • Balance (minimal value: 50 mg; the precision requirement: 50 ± ≤0.001 mg) and the
98 appropriate mass standards F1, e.g. 1 - 500 mg ± ≤0.02 mg; 1 g ± ≤0.01 mg
- 99 • Clean glass vials for the preparation of stock solutions (e.g., gas chromatography vials or
100 Sarstedt #86.1509)
- 101 • 96-well microplates (e.g., Nunc #167008)
- 102 • Freezer below -16 °C
- 103 • Refrigerator at 2-10 °C
- 104 • CO₂ humidified incubator at 37°C +/- 2 °C, 5% CO₂ +/- 0.5%
- 105 • 37 °C water bath

- 106 • PipetteAid
- 107 • Pipettors (p1000, p200, p20) or Micropipettes (1000-100; 100-10; 10-0.5)
- 108 • Centrifuge
- 109 • Vacuum aspirator
- 110 • Microplate Reader; for measuring absorbance
- 111 • Laminar Flow Hood
- 112 • Culture Flasks (T75, T25; e.g. Nunc, Falcon)
- 113 • Serological pipettes (1 mL, 5 mL, 10 mL, and 25 mL)
- 114 • Sterile, filter pipette tips (10-20 µL, 200 µL, and 1000 µL)
- 115 • Sterile, pipette tips without filter (e.g. 10-20 µL, 200 µL, and 1000 µL)
- 116 • Conical tubes (e.g. 5 mL, 15 mL, and 50 mL)
- 117 • Polypropylene Cluster Tubes (1.1 mL e.g. Corning #MTS-11-8-C)
- 118 • Syringe filters (0.22µm)

119 **1.6 Reagents and Chemicals**

- 120 • Foetal calf/bovine serum (e.g. Gibco #10270106)
- 121 • Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12; 1:1 mixture with phenol red and HEPES (15 mM) (e.g. Sigma-Aldrich #D6421)
- 122 • DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES and w/o phenol red (e.g. Sigma-Aldrich #D6434)
- 123 • HEPES 1M solution (e.g. Sigma-Aldrich #H0887)
- 124 • Penicillin-Streptomycin Solution (e.g. Sigma-Aldrich #P0781)
- 125 • Bovine insulin (e.g. Sigma-Aldrich #I0516)
- 126 • Ethanolamine (e.g. Sigma-Aldrich #E0135)
- 127 • Sodium selenite (e.g. Sigma-Aldrich #S5261)
- 128 • Human apotransferrin (e.g. Sigma-Aldrich #T2036)
- 129 • Bovine serum albumin (e.g. Sigma-Aldrich #A9418)
- 130 • DMSO 99.9% purity (e.g. Sigma-Aldrich #D8418)
- 131 • Ethyl alcohol 99.8% pure p.a. (e.g. POCH #396480111)
- 132 • Deionised water (dH₂O)

136 **1.7 Reference items, positive and negative control Item(s)**

	Agonism	Antagonism
Reference item	3,3'-5-triiodothyronine (T3) [CAS 6893-02-3], ≥ 95% purity	5,5-Diphenylhydantoin (DPH, or Phenytoin) [CAS 57-41-0]
Positive control item	T4, L-Thyroxine (3,3',5,5''-Tetraiodo-L-thyronine) [CAS 51-48-9], ≥ 98%	
Negative control item	Mefenamic acid (MfA) [CAS 61-68-7]	
Cytotoxic positive control item	Sodium Dodecyl Sulfate (SDS) [CAS 151-21-3]	

137

138 **1.8 Preparations of media and reagents**

139 **1.8.1 Cell culture medium (cDMEM/F12; completed cDMEM/F12)**

- 140 • Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12; 1:1
141 mixture; with phenol red and 15 mM HEPES) supplemented with 10% Foetal calf serum, 2.5
142 mM L-Glutamine and additional amount of HEPES to obtain 25 mM.

143 To prepare cDMEM/F-12 the following volume of supplements should be added:

	Final concentration	Volume [mL]
DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES (e.g. Sigma-Aldrich #D6421)*		500
Heat inactivated FCS/FBS**	10%	57.6
200 mM L-Glutamine	2.5 mM	7.2
1 M HEPES	up to 25 mM	5.8
Penicillin-Streptomycin solution (100x)	1%	5.8
Total		576.4

144 *if different DMEM/F-12 is used the final concentration of HEPES and L-Glutamine should be
145 adjusted to 2.5 mM and 25 mM, respectively.

146 ** if heat-inactivated FCS/FBS was purchased, it should be defrost at 2-10°C e.g. during the night and
147 stored in aliquots at below -16°C.

148 If non-inactivated FCS/FBS was purchased, it should be heat-inactivated in water bath at 56-57°C for
149 30-35 min and filtered (0.2 µm) before aliquoting.

150 **1.8.2 PCM medium**

151 PCM medium is a Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12;
152 1:1 mixture) without phenol red supplemented with 10 µg/mL bovine insulin, 10 µM ethanolamine,
153 10 ng/mL sodium selenite, 10 µg/mL human apotransferrin, 500 µg/mL bovine serum albumin.
154 Neither T3 nor T4 are added to the PCM medium.

155 To prepare PCM medium the following volume of supplements should be added:

	Final concentration	Volume [mL]
DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES (e.g. Sigma-Aldrich # D6434)*		500
200 mM L-Glutamine	2.5 mM	6.5
1 M HEPES	up to 25 mM	5.2
Penicillin-Streptomycin solution (100x)	1%	5.2
10 mg/mL bovine insulin	10 µg/ml	0.523
50 mM ethanolamine (Section 1.8.5-1)	10 µM	0.105
50 µg/mL sodium selenite (Section 1.8.5-2)	10 ng/mL	0.105
2 mg/mL human apotransferrin (Section 1.8.5-3)	10 µg/mL	2.6
100 mg/mL bovine serum albumin (Section 1.8.5-4)	500 µg/mL	2.6
Total		523.0

156 *if different DMEM/F-12 is used the final concentration of HEPES and L-Glutamine should be
157 adjusted to 2.5 mM and 25 mM, respectively.

158 **1.8.3 PCM medium +4x EC₅₀ concentration of T3 (PCM-T3 medium)**

159 Once for each working cell bank and for each new batch of T3 the laboratory should verify the EC50
160 value for reference item T3 as described in Section 2. If the EC50 value of T3 meets the acceptance
161 criteria (Section 2.4) then the T3 concentration of 0.1 nM should be used as the EC50 value, i.e. PCM
162 medium with 0.4 nM T3 will be considered as PCM-T3 medium (4x EC50).

163 For each test item 2.5 mL of PCM-T3 medium is needed. To this end, prepare 3-step dilution of T3
164 stock solution (1 mM) in PCM medium to get 0.4 nM of T3 (1 mM→(100x)→10 μM→(100x)→100
165 nM→(250x)→ 0.4 nM). Since the solvent concentration in PCM-T3 is lower than 0.001% (exactly
166 0.00004%) and the maximum accepted concentration for DMSO is 0.5%, the effect of the solvent is
167 omitted (the final concentration of DMSO in test or control items during the T-screen test is
168 0.40004%). PCM-T3 medium will be used for antagonist experiments – Section 3.2-2c and Section
169 3.3.3-2.

170 **1.8.4 Stock and working solutions of chemicals**

171 To prepare the **exposure concentration** of every chemical given below, firstly the **concentrated**
172 **solution** is prepared that is used to prepare **the working solution** (the double desired exposure
173 concentration) according to Section 3.2-2b and/or Section 3.2-2c.

174 1.8.4-1 **3,3'-5-triiodothyronine (T3) [1mM stock solution] and [4 nM working solution]**

175 Dissolve T3 in DMSO to produce a 1 mM stock solution. Use fresh or store in aliquots below
176 -16 °C.

177 T3 is used in the T-Screen assay in the full dose range (the first Agonism plate) and in the
178 highest test concentration, i.e. 2 nM – to prepare the *concentrated solution* of T3 [1 μM],
179 dilute T3 stock solution in the solvent/DMSO as follows: 1 mM → (50x) → 20 μM --(20x)→
180 1 μM. Then, use the *concentrated solution* of T3 to prepare *working solution* according to
181 Section 3.2-2b (only in Agonism plate). T3 is added at the final concentration EC50 into
182 every test well on Antagonism plates except UC (PCM) and the preparation of sample with
183 addition of T3 EC50 is described in Section 1.8.3.

184 1.8.4-2 **3,3',5,5'-tetraiodothyroxine (T4) [2 mM stock solution] and [20 nM working solution]**

185 Dissolve T4 in DMSO to produce a 2 mM stock solution. Use fresh or store in aliquots below
186 -16 °C.

187 The exposure concentration of T4 in the T-Screen assay is 10 nM – to prepare the
188 *concentrated solution* of T4 [5 μM], dilute T4 stock solution in the solvent/DMSO as
189 follows: 2 mM → (400x) → 5 μM. Then, use *the concentrated solution* of T4 to prepare
190 *working solution* according to Section 3.2-2b, for Agonism experiment.

191 1.8.4-3 **5,5-Diphenylhydantoin (Phenytoin; DPH) [50 mM stock solution] and [100 μM working solution]**

192 Dissolve DPH in DMSO to produce a 1 mM stock solution. Use fresh or store in aliquots
193 below -16 °C.

194 The exposure concentration of DPH in the T-Screen assay is 50 μM – to prepare the
195 *concentrated solution* of DPH [25 mM], dilute DPH stock solution in the solvent/DMSO as
196 follows: 50 mM → (2x) → 25 mM . Then, use *the concentrated solution* of DPH to prepare
197 *working solution* according to Section 3.2-2c, for Antagonism experiment.

198 1.8.4-4 **Mefenamic acid (MfA) [100 mM stock solution] and [200 nM working solution]**

200 Dissolve MfA in DMSO to produce a 100 mM stock solution. Use fresh or store in aliquots
201 below -16°C.

202 The exposure concentration of MfA in the T-Screen assay is 100 nM – to prepare the
203 *concentrated solution* of MfA [50 µM], dilute MfA stock solution in the solvent/DMSO as
204 follows: 100 mM → (100x) → 1 mM --(20x)→ 50 µM. Then, use the *concentrated solution*
205 of MfA to prepare *working solution* according to Section 3.2-2b and Section 3.2-2c, for
206 Agonism and Antagonism experiments, respectively.

207 1.8.4-5 **Sodium Dodecyl Sulfate (SDS) [100 mM stock solution] and [200 µM working solution]**

208 Prepare the 100 mM stock solution of SDS in DMSO. Use fresh or store in aliquots at room
209 temperature.

210 The exposure concentration of SDS in the T-Screen assay is 100 µM – to prepare the
211 *concentrated solution* of SDS [50 mM], dilute SDS stock solution in the solvent/DMSO as
212 follows: 100 mM → (2x)→ 50 mM. Then, use the *concentrated solution* of SDS to prepare
213 working solution according to Section 3.2-2b (SDS is used only on Agonism plate).

214 **1.8.5 Stock solutions of reagents**

215 1.8.5-1 **Ethanolamine [50 mM]**

216 Prepare a 50 mM stock solution of ethanolamine in dH₂O by diluting a pure ethanolamine
217 (16.6 M) 332x, i.e. 10 µL of ethanolamine (16.6 M) added to 3.310 mL of dH₂O, filter (0.22
218 µm pore size). Use fresh or store at room temperature in closed glass containers.

219 1.8.5-2 **Sodium selenite [50 µg/mL]**

220 Prepare a stock solution of sodium selenite (50 µg/mL) in dH₂O, filter (0.22 µm pore size)
221 and store in aliquots below -16 °C.

222 1.8.5-3 **Human apotransferrin [2 mg/mL]**

223 Prepare a 2 mg/mL stock solution of human apotransferrin in dH₂O, filter (0.22 µm pore
224 size), and store in aliquots below -16 °C. Stock solutions of apo-transferrin stored at 2-8°C
225 are stable for 5–10 days.

226 1.8.5-4 **Bovine serum albumin [100 mg/mL]**

227 Prepare a 100 mg/mL stock solution of bovine serum albumin in dH₂O, filtered (0.22 µm
228 pore size) and store in aliquots at 2-8°C

229 **2. Verification of the EC50 value of T3**

230 Once for each working cell bank and for each new batch of T3 the laboratory should verify the EC50
231 value for reference item T3 as described below.

232 **2.1 Seeding the GH3 cells onto 96-well plate**

233 Perform the cell seeding as described for the T-screen (Section 3.3.1). Seed cells onto a 96-well
234 microplate in three replicates (e.g. rows B-D; Figure 1).

235 **2.2 Exposure of the cells to the range of T3 concentrations**

236 2.2 - 1. Prepare T3 at the range of double desired concentrations by diluting the concentration
237 of 4 nM using dilution factor (DF) = 3. The final 7 concentration of T3 should be: 2;
238 0.667; 0.222; 0.074; 0.025; 0.008 and 0.003 nM

- 239 2.2 - 2. Add 100 µL T3 in PCM medium at the appropriate concentration into the 100 µL already
 240 present (with the cells) in the respective well(s) (to achieve a 1:1 dilution). Test T3 and
 241 control items in triplicate (Figure 1).
 242 2.2 - 3. Incubate the plate for 96 ± 1 h at 37 ± 2 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid
 243 atmosphere.
 244

245 **Figure 1. Scheme of the plate layout for verification of the EC50 value of T3**

	1	2	3	4	5	6	7	8	9	10	11	12
A	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z
B	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
C	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
D	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
E	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z
F												
G												
H												

246 **SC** (wells 2B-D and 11B-D) - solvent control

247 **PCM** (wells 10B-D) - test system control

248 **T3 C** (wells 3B-9D) - the range of T3 concentrations (0.003 - 2 nM); where C1 is the lowest T3
 249 concentration

250 **Z** – external wells (medium or PBS)

251 **Blank** - PCM medium without cells (AlamarBlue is added in the proliferation assay); used for
 252 calculations both % AlamarBlue reduction (%AR) and % Dye reduction (%DR) according to
 253 formulas given in SOP *Determination of cell proliferation in T-screen assay* in Section 2.2.1 or
 254 Section 2.2.2, respectively

255 **Blank w/o AB** - PCM medium without cells (AlamarBlue is not added in the proliferation assay);
 256 used for calculations only %DR according to formula given in SOP *Determination of cell*
 257 *proliferation in T-screen assay* in Section 2.2.2.

258 2.3 Measurement of cell proliferation

259 2.3 - 1. Perform the cell proliferation assay according to Section 2: “AlamarBlue assay”
 260 described in SOP “*Determination of cell proliferation in T-screen assay*”.

261 2.3 - 2. Based on received results (optical density/absorbance), calculate **% AlamarBlue**
 262 **reduction (%AR)** or **% Dye reduction (%DR)** for each triplicate sample using formulas
 263 given in SOP *Determination of cell proliferation in T-screen assay* in Section 2.2.1 or
 264 Section 2.2.2, respectively.

265 2.3 - 3. Calculate the increase of *Cell proliferation* that is expressed as **the relative proliferative**
 266 **effect (RPE)** according to the formula given in Section 3.4.1-2.

267 2.3 - 4. Calculate the **EC50** of T3 according to the formula given in Section 3.4.1-3.

268 2.4 Acceptance criteria

269 Based on the obtained results the mean EC50 value of T3 should be calculated. The mean EC50
 270 value should be $-10 \pm 0.4 \log_{10}(\text{Molar})$ units (in the range from -10.4 to $-9.6 \log_{10}(\text{Molar})$ units).


271 If acceptance criteria for the mean EC50 value of T3 are met, 0.1 nM T3 should be used as the
 272 EC50 value of T3 in the T-screen test.

273 Every result, i.e. the mean EC50 of T3, should be registered to create the historical data EC50 value
 274 of T3. Also, the mean EC50 of T3 calculated from the results obtained during every T-screen test (the
 275 first agonism plate, Section 3.3.3) should be included to the historical data.

276 **3. T-Screen experimental procedure**

277 Figure 2 summarises the main steps required to run this method.

278 **Figure 2.** The T-screen assay workflow with absorbance of AlamarBlue dye detection.



279

280 **3.1 Preparation of Test Item**

281 The test item's solubility in the solvent and medium will need to be assessed. Both stock solutions
 282 and working solutions should be assessed.

283 **3.1.1 Stock solutions**

284 *It is recommended to use the same solvent for all reference and control items and to use that solvent*
 285 *for the test items, to ensure all samples are tested under the same conditions. In case another solvent*
 286 *is used for the test item, the impact on the test system and the results must be assessed.*

287 3.1.1-1 Prepare a 100 mM or 10 mg/mL stock solution (may be adapted on basis of data if too low
 288 or high) of test item by weighing a nominal amount into a clear glass vial and add the
 289 appropriate volume of solvent (use DMSO as a default solvent and if not possible to dissolve
 290 test item in DMSO then other solvents, e.g. water, PCM medium or ethanol should be used.
 291 For visual inspection, it is important to have at least 0.5 ml in the vial to be able to observe
 292 well.

293 3.1.1-2 Visually check if dissolved. If not dissolved, heat (37°C; up to 30 min) may be applied to aid
 294 solubility.

295 3.1.1-3 If not possible to solubilise, prepare a lower (e.g. 1:2) stock solution (or dilute existing
 296 stock) and again check solubility. Change solvent if needed.

297 3.1.1-4 Continue until the stock solution is soluble.

298 **3.1.2 Work solutions to determine solubility**

299 3.1.2-1 Dilute the stock solution of test item in solvent prepared in Section 3.1.1. Prepare the
 300 range of 7 concentrations (*the concentrated solutions*) using a dilution factor 10 (10 µL of
 301 stock solution of the appropriate *concentrated solution* to 90 µL solvent in a 96-well plate).

302 For test item dissolved in DMSO or ethanol prepare the 500x *concentrated solutions*,
303 whereas for test item dissolved in water or PBS, prepare 100x *concentrated solutions*.

304 3.1.2-2 Prepare the 2x concentrated solutions (*working solutions*) in PCM medium by transferring:
305 2 µL of the appropriate *concentrated solutions* to 498 µL of PCM in a 24-well plate or 1.1
306 mL cluster tubes.

307 Solvent concentration in the *working solutions* of test item should not exceed 0.4% for
308 DMSO and ethanol or 10% for water or PBS (the final concentration 0.2% or 5%,
309 respectively for DMSO/ethanol or water/PBS).

310 3.1.2-3 Visually check if the working solutions are dissolved.

311 3.1.2-4 Identify the highest soluble concentration.

312 **3.2 Pre-screen experiment**

313 During the pre-screen experiment it is assessed if the test item is:

- 314 • A Thyroid Hormone Agonist
- 315 • A Thyroid Hormone Antagonist
- 316 • Cytotoxic

317
318 and at which dose-range it increases or decreases cell proliferation (Range finder).

319
320 **Note:** GH3 cells are not able to divide properly in PCM without T3 but keep basal or low
321 activity. Concentrations of test items that lower the cellular activity of GH3 cells cultured in
322 PCM medium without T3 (as determined with the cell proliferation assay) are considered to be
323 cytotoxic.

324 **3.2.1 Seeding the GH3 cells onto 96-well plate**

325 Perform the cell seeding as described for the T-screen (Section 3.3.1).

- 326 - Seed cells into all internal wells, i.e. B2-G11.
- 327 - Add 100 µL of PCM medium to each well in columns 1 and 12, for *Blank* and *Blank w/o AB*
328 (Figure 4) that are used for calculations in the proliferation assay.
- 329 - Add 100 µL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
330 proper humidity for the cells.

331 **3.2.2 Preparation of test, reference and control item solutions**

332 On the day of treatment, prepare the range of 7 *working solutions*, both for **Agonist experiments**
333 and **Antagonist assessment**, as follows:

334 3.2.2-1 Stage 1: Prepare stock solution of test item at the concentration and solvent determined
335 above (Section 3.1.1). Prepare a dilution series of test item stock in the solvent (the range of
336 7 *concentrated solutions*) as described in Section 3.1.2-1 (Figure 3; Stage 1).

337 In the first experiment a dose range finding is performed to capture the whole dose response range
338 for the next experiments. The highest concentration to be tested is the highest non cytotoxic soluble
339 concentration and dilution factor 10 is applied for the endpoint measurement. In the next
340 experiment a dose range includes 7 concentrations where the highest concentration shows the max
341 effect (the induction of proliferation in Agonist experiment) and the lowest concentration shows no

342 effect (effect comparable to solvent control in Agonist experiment). Dilution factor is adjusted to
343 cover 7 concentrations and each test item have its own dilution factor (e.g. DF 1.5, 2, 3, 4 or 5).

344 Use the same range of the 7 *concentrated solutions* to prepare *working solutions*, both for **Agonist**
345 (Stage 2a) and **Antagonist assessment** (Stage 2b).

346
347 3.2.2-2 For **Agonist assessment (Treatment without T3)** (Figure 3; Stage 2a)

348 Stage 2a: Prepare *working solutions* as described in Section 3.1.2-2.

349 In the first experiment/the dose range finding an interference of test item with the
350 assay/AlamarBlue reagent should be assessed. To this end, prepare *working solutions* as described in
351 Section 3.1.2-2 using double volumes to have enough solutions both for test plate and the additional
352 plate for testing interference of AlamarBlue with test item.

353

354 3.2.2-3 For **Antagonist assessment – Treatment with EC₅₀ of T3** (Figure 3; Stage 2b)

355 Stage 2b: Prepare *working solutions* containing 2xEC₅₀ of T3 (dilute *concentrated solutions* prepared
356 in Stage 1 in the appropriate mixture of PCM medium and PCM-T3 medium), as follows:

- 357 - Add 248 µL PCM medium into 1.1 mL-cluster tubes, then
- 358 - Add 250 µL PCM-T3 medium (4x EC₅₀ of T3; prepared as described in Section 1.8.3) into 1.1
359 mL-cluster tubes, then
- 360 - Add 2 µL of the appropriate *concentrated solution* to the mixture of PCM medium and PCM-
361 T3 medium already present into 1.1 mL-cluster tubes. Solvent concentration (both for test
362 item and T3) should not exceed 0.5%.

363 3.2.3 Exposure of the cells

364 Perform the exposure as follows:

365 3.2.3-1 Upper part of the plate/B2:D11 (agonism assessment):

- 366 - Add 100 µL *working solutions (the double desired exposure concentrations)* of test item
367 (TI) or solvent control (SC) prepared in PCM medium to the 100 µL already present (with
368 the cells) in the respective well(s) (to achieve a 1:1 dilution). Test all samples in triplicate.
- 369 - Add 100 µL PCM medium only to three wells designed UC (PCM) (untreated control/cell
370 system control; Figure 4)).
- 371 - Add the same samples in the same order and volume (100 µL) into the additional plate
372 prepared for testing interference of AlamarBlue with test item (Section 3.2-3e).

373 3.2.3-2 Lower part of the plate/E2:G11 (antagonism assessment):

- 374 - Add 100 µL test item or solvent control (S/T3) *working solutions (the double desired*
375 *exposure concentrations)* prepared in PCM medium with addition of 2x EC₅₀ value of T3
376 to the 100 µL already present (with the cells) in the respective well(s) (to achieve a 1:1
377 dilution). Test all samples in triplicate.
- 378 - Add 100 µL PCM medium with addition of 2x EC₅₀ value of T3 to three wells designed
379 (EC₅₀ T3; Figure 4).
- 380 - Add 100 µL 2x Ref(T3)C1 = the highest concentration in the range described in Section
381 3.1.2.
- 382 - Add the same samples in the same order and volume (100 µL) into the additional plate
383 prepared for testing interference of AlamarBlue with test item (Section 3.2.3-5).

384 3.2.3-3 Add 100 µL of PCM medium to the 100 µL already present (without the cells) in each well
385 in columns 1 and 12 (Blank and Blank w/o AB; Figure 4) that are used for calculations in the
386 proliferation assay.

387 3.2.3-4 Add 100 µL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
388 proper humidity for the cells.

389 3.2.3-5 Prepare the additional plate for testing interference of AlamarBlue reagent with test items
390 (samples incubated in PCM medium without cells) – add 100 µL of PCM medium into every
391 well needed (triplicates/sample). Then, add all prepared samples (100 µL) in the same
392 order as indicated in Section 3.2.3-1 and Section 3.2.3-2. Also, add extra PCM medium as
393 given in Section 3.2.3-3 and Section 3.2.3-4.

394 3.2.3-6 Incubate plates for 96 ± 1 h at 37 ± 2 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid atmosphere.
395 *Plate layout for Pre-screen experiment is presented in Figure 4. It is advised to assess agonism and*
396 *antagonism effect on the same plate (both in triplicates). For antagonism all treatments are*
397 *performed in presence of T3.*

398 **3.2.4 Measurement of cell proliferation and data analysis**

399 3.2.4-1 Perform the AlamarBlue assay according to Section 2 in SOP "Determination of cell
400 proliferation in T-screen assay" to assess cell proliferation or cytotoxic effects.

401 3.2.4-2 Based on received results (optical density/absorbance), calculate % AlamarBlue reduction
402 (%AR) or % **Dye reduction (%DR)** for each triplicate sample using formulas given in SOP
403 *Determination of cell proliferation in T-screen assay* in Section 2.2.1 or Section 2.2.2,
404 respectively.

405 3.2.4-3 Calculate the **Cell proliferation (%CP)** for each test item (the concentration range) based on
406 values %AR and %DR calculated above (Section 3.2.4-2) according to the following formulas,
407 respectively for calculated %DR or %AR:

$$\%CP = (\%DR_{x/T3} - \%DR_x) \times 100$$

408 **OR***

$$\%CP = (\%AR_{x/T3} - \%AR_x) \times 100$$

409 where:

410 x - the effect of TI, PC, NC or REF T3 at the concentration analysed; tested without EC50 of T3

411 x/T3 - the effect of TI, PC, NC or REF T3 at the concentration analysed; tested in the presence of
412 EC50 of T3

413 * - %AR or %DR is used depending on the method choose for calculation results of the
414 proliferation assay (SOP *Determination of cell proliferation in T-screen assay* in Section
415 2.2.1 or Section 2.2.2, respectively)

416 3.2.4-4 Calculate the increase of cell proliferation that is expressed as **the relative proliferative**
417 **effect (RPE)** for "EC50 T3" according to Section 3.4.1-2.

418

419 **3.2.5 Acceptance criteria for pre-screen assay**

420 3.2.5-1 %DR or %AR for UC(PCM) should not be more than 15% different from TI SC

421 3.2.5-2 RPE for EC50 T3 should be in the range of 30-70%

422 **3.2.6 Identification of the range concentrations of TI for the next experiments**

423 Identify the range concentrations of TI for the next experiments that captures the whole dose
424 response range (from no effect to the highest agonistic/antagonistic effect). The range of
425 concentration should have not more than two concentrations given the maximal response and not
426 more than two concentrations given the minimal response (near background/ not more than 5% of
427 SC).

428 3.2.6-1 If a strong cytotoxic effect is detected, i.e. three or more concentrations are found to be
429 cytotoxic, repeat the test using the changed range of concentrations (start from the highest
430 non-cytotoxic concentration and adjust DF to get 7 concentrations) to better identify range
431 of test concentrations of test item. The highest concentration for the T-screen assay should
432 be the highest noncytotoxic concentration of test item.

433 3.2.6-2 If no cytotoxic effect is observed, but the agonistic or antagonistic effect is not observed
434 choose for the next experiment the same range of concentration and DF 10

435 3.2.6-3 If no cytotoxic effect is observed, but any agonistic or antagonistic effect is observed,
436 choose DF to obtain the whole dose response (e.g. DF 1.5, 2, 3, 4 or 5) (if any).


437

438

CONFIDENTIAL

439 **Figure 3.** A schema of test item preparation for pre-screen experiment and the T-Screen test.


440 **Stage 1. Dilution of a stock solution**



The solvent concentration in the diluted stock solutions/*concentrate solutions* is considered as 100%.

Stage 2a. Preparation of the *working solutions* (double desired exposure concentrations) for AGONIST experiments.


Cells should be exposed to a test item from minimal to maximal concentration.



The solvent concentration in the *working solutions* is considered as 0.4 %.

Stage 2b. Preparation of the *working solutions* (double desired exposure concentrations) for ANTAGONIST experiments.

Cells should be exposed to a test item from minimal to maximal concentration.



Since the solvent concentration in PCM-T3 medium is only 0.00004%, the effect of the solvent is omitted. The total solvent concentration in the *working solutions* can be considered as 0.4% - regardless of whether solvent of test item is the same as solvent of T3 (i.e. DMSO) or different e.g. ethanol.

441 **Figure 4.** Plate layout for pre-screen experiment

	1	2	3	4	5	6	7	8	9	10	11	12	
A	*	*	*	*	*	*	*	*	*	*	*	*	
agonism	B	Blank	SC	TI C7	TI C6	TI C5	TI C4	TI C3	TI C2	TI C1	UC (PCM)	SC	Blank w/o AB
	C	Blank	SC	TI C7	TI C6	TI C5	TI C4	TI C3	TI C2	TI C1	UC (PCM)	SC	Blank w/o AB
	D	Blank	SC	TI C7	TI C6	TI C5	TI C4	TI C3	TI C2	TI C1	UC (PCM)	SC	Blank w/o AB
antagonism	E	Blank	S/T3	TI C7/ T3	TI C6/ T3	TI C5/ T3	TI C4/ T3	TI C3/ T3	TI C2/ T3	TI C1/ T3	EC50 T3	Ref(T3)C1	Blank w/o AB
	G	Blank	S/T3	TI C7/ T3	TI C6/ T3	TI C5/ T3	TI C4/ T3	TI C3/ T3	TI C2/ T3	TI C1/ T3	EC50 T3	Ref(T3)C1	Blank w/o AB
	G	Blank	S/T3	TI C7/ T3	TI C6/ T3	TI C5/ T3	TI C4/ T3	TI C3/ T3	TI C2/ T3	TI C1/ T3	EC50 T3	Ref(T3)C1	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*	

442

UC (PCM)	Untreated Control; Tested as triplicate sample on each plate, PCM Medium Only
EC50 T3	PCM medium + EC50 T3; Tested as triplicate sample on each plate
SC	PCM medium + solvent; Tested as double triplicate sample on each side of the plate
S/T3	PCM medium + solvent + EC50 T3; Tested as double triplicate sample on each side of the plate
Ref(T3)C1	the max. concentration of T3 used in the study/the concentration no. 1; 2 nM); it will be used for calculation of RPE for EC50
*	Outer wells with PBS or PCM medium only; Outer wells may only be used only when the plate has the additional collar/space to fill it with liquid to ensure the cells proper humidity
TI C [7-1]	Test item in the range of concentration (C) from C7/Cmin to C1/Cmax (C7 is the lowest concentration tested) TI tested without addition of T3 (agonism experiment)
TI C [7-1] / T3	Test item in the range of concentration (C) from C7/Cmin to C1/Cmax (C7 is the lowest concentration tested); TI tested in the presence of EC50 T3 (antagonism experiment)
Blank	PCM medium without cells (AlamarBlue is added in the proliferation assay); used for calculations both % AR and % DR according to formulas given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.2.1 or Section 2.2.2, respectively
Blank w/o AB	PCM medium without cells (AlamarBlue is not added in the proliferation assay); used for calculations only %DR according to formula given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.2.2

443 **3.3 T-Screen assay**

444 IMPORTANT! Agonistic and antagonistic potential of test item should be assessed simultaneously but
445 on separate plates. Complete Reference item standard curves (for Ref(T3) and Ref(DPH) in the
446 agonism and antagonism plates, respectively) should be included in each assay.

447 **3.3.1 Seeding the GH3 cells onto 96-well plate**

448 Number of plates depends on number of test items. Because agonistic and antagonistic potential of
449 TI should be assessed on separate plates, to test one TI two plates are needed. Because complete
450 Reference item standard curves should be included in each assay thus the first set of two plates is
451 used to assess one TI together with REF (T3) or REF (DPH) on the agonism or antagonism plates,
452 respectively. The second set of two plates is used to assess the next two TI (Figure 5).

453 3.3.1-1 Forty-eight hours prior to plating the cells onto 96-well microplates for the experiment,
454 change the standard culture medium to PCM medium (Section 1.8.2).

455 3.3.1-2 Release cells as described in SOP "Handling and Maintenance of GH3 cell line", Section
456 2.3.4 "Detachment of GH3 cells". Pipet cell suspension very carefully several times (cells
457 easily detach but are sensitive to shaking as well as the effects of trypsin).

458 3.3.1-3 Determine density of cell suspension taking into account viable cells as described in SOP
459 "Handling and Maintenance of GH3 cell line", Section 2.3.6.

460 3.3.1-4 Calculate cell viability as described in SOP "Handling and Maintenance of GH3 cell line",
461 Section 2.3.6.

462 3.3.1-5 Only if viability of GH3 is more than 90%, dilute cell suspension in PCM medium and seed
463 cells onto a 96-well microplate at a density of 2500 viable cells/well in 100 µL PCM
464 medium.

465 3.3.1-6 Pre-incubate for 24 hours +/- 2 hours at 37 ± 2 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid
466 atmosphere to allow cells to attach to bottom of wells before the treatment.

467 3.3.1-7 Before the treatment, based on microscopic observation, the following acceptance
468 criterion should be met (otherwise, the plate is rejected): not more than 50% of cells in the
469 well are floating

470 **3.3.2 Preparation of test, reference and control item solutions**

471 Prepare test item *working solutions* as for the pre-screen experiment, using the dilution factor (DF)
472 identified in the range finder pre-screen experiment that will capture the whole dose response. If no
473 agonistic/antagonistic effect is observed the same range of concentrations (and the same DF) should
474 be used in T-screen assay. If any agonistic/antagonistic effect is observed adjust DF to obtain the
475 whole dose response (e.g. DF 2, 3, 4 or 5) (if any) and then repeat the experiment.

476 Reference and control items are used in T-Screen assay using the fixed concentration – the
477 preparation of *working solutions* of reference and control items (*the double desired exposure*
478 *concentrations*) are described in Section 1.8.4.

479 **3.3.3 Exposure of the cells**

480 *Plate layout for the 8a T-Screen assay is presented in Figure 5.*

481 For every test item the set of two plates is prepared, as follows:

482 3.3.3-1 **Agonist plates:**

- 483 - Add 100 μL *working solutions (the double desired exposure concentrations)* of test item (T1) or
484 the appropriate controls prepared in PCM medium to the 100 μL already present (with the
485 cells) in the respective well(s) (to achieve a 1:1 dilution). Test all samples in triplicate with
486 exception of solvent control (SC) that is tested in 6-12 repetitions spread out each side of the
487 plate.
- 488 - Add 100 μL 2x Ref(T3)C1 = the highest concentration in the range described in Section 3.1.2.
- 489 - Add 100 μL PCM medium to three wells designed UC (PCM) (untreated control/cell system
490 control; Figure 5).
- 491 - Add 100 μL of PCM medium to the 100 μL already present (without the cells) in each well in
492 columns 1 and 12 (*Blank* and *Blank w/o AB*; Figure 5) that are used in the proliferation assay
493 for calculations.
- 494 - Add 100 μL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
495 proper humidity for the cells.
- 496 - Incubate plates for 96 ± 1 h at 37 ± 2 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid atmosphere.

497

498 3.3.3-2 **Antagonism plates:**

- 499 - Add 100 μL *working solutions (the double desired exposure concentrations)* of test item
500 (T1/T3) or the appropriate controls prepared in PCM medium with addition of 2x EC50 value
501 of T3 to the 100 μL already present (with the cells) in the respective well(s) (to achieve a 1:1
502 dilution). Test all samples in triplicate with exception of solvent control (S/T3) that is tested
503 in 6-9 repetitions spread out each side of the plate.
- 504 - Add 100 μL 2x Ref(T3)C1 = the highest concentration in the range described in Section 3.1.2.
- 505 - Add 100 μL PCM medium (without EC50 T3!) to three wells designed UC (PCM) (untreated
506 control/cell system control; Figure 5).
- 507 - Add 100 μL of PCM medium to the 100 μL already present (without the cells) in each well in
508 columns 1 and 12 (*Blank* and *Blank w/o AB*; Figure 5) that are used in the proliferation assay
509 for calculations.
- 510 - Add 100 μL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
511 proper humidity for the cells.
- 512 - Incubate plates for 96 ± 1 h at 37 ± 2 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid atmosphere.

513 **Figure 5.** Plate layout for the 8a T-Screen assay.

514 (A) the first set of plates consists of **the agonism plate 1** that is used to assess the range of Ref (T3) concentrations (upper part of the plate) and test item 1
 515 (lower part of the plate) together with all appropriate controls and **the antagonism plate 1** that is used to assess the range of Ref (DPH) concentrations
 516 (upper part of the plate) and test item 1 (lower part of the plate) in the presence of EC50 T3 together with all appropriate controls
 517

REF (T3)	Reference item T3 for AGONISM; Tested as full dose response curve once per series of plates. Each plate should contain the triplicate samples used in data analysis on the plate
REF (DPH)	Reference item Amiodarone for ANTAGONISM; Tested as full dose response curve once per series of plates. Each plate should contain the triplicate samples used in data analysis on the plate
.../T3	Sample tested in the presence of EC50 T3
PC (A)	Positive control for AGONISM; C max; Tested as triplicate sample once per series of plates
PC (ANT)	Positive control for ANTAGONISM; C max; Tested as triplicate sample once per series of plates
NC	Negative Control; C max; Tested as triplicate sample once per series of plates; the same chemical for AGONISM and ATAGONISM
UC (PCM)	Untreated Control; Tested as triplicate sample on each plate, PCM Medium Only,
SC	PCM medium + solvent; Tested as triplicate sample on each plate. It is recommended to use the same solvent for all reference and control items and to use that solvent for the test items, to ensure all samples are tested under the same conditions. In case another solvent is used for the test item, the impact on the test system and the results must be assessed.
SDS	Cytotox control; Tested as triplicate sample once per series of plates; only tested on AGONISM plate
*	Outer wells with PBS or PCM medium only; Outer wells may only be used only when the plate has the additional collar/space to fill it with liquid to ensure the cells proper humidity.
TI	Test item; Tested as full dose response curve
C	Concentrations of TI or REF (from C7/Cmin to C1/Cmax; C7 is the lowest concentration tested)
Blank	PCM medium without cells (AlamarBlue is added in the proliferation assay); used for calculations both % AR and % DR according to formulas given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.2.1 or Section 2.2.2, respectively
Blank w/o AB	PCM medium without cells (AlamarBlue is not added in the proliferation assay); used for calculations only %DR according to formula given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.2.2

518

519 **Figure 5. continued**

520

521 **(A) the first set of plates:**

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
C	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
D	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
E	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
G	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
G	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

522

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	S/T3	Ref(DPH)C7/T	Ref(DPH)C6/T	Ref(DPH)C5/T	Ref(DPH)C4/T	Ref(DPH)C3/T	Ref(DPH)C2/T	Ref(DPH)C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
C	Blank	S/T3	Ref(DPH)C7/T	Ref(DPH)C6/T	Ref(DPH)C5/T	Ref(DPH)C4/T	Ref(DPH)C3/T	Ref(DPH)C2/T	Ref(DPH)C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
D	Blank	S/T3	Ref(DPH)C7/T	Ref(DPH)C6/T	Ref(DPH)C5/T	Ref(DPH)C4/T	Ref(DPH)C3/T	Ref(DPH)C2/T	Ref(DPH)C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
E	Blank	NC/T3	TI 1 C7/T3	TI 1 C6/T3	TI 1 C5/T3	TI 1 C4/T3	TI 1 C3/T3	TI 1 C2/T3	TI 1 C1/T3	PC(ANT)/T3	S/T3	Blank w/o AB
G	Blank	NC/T3	TI 1 C7/T3	TI 1 C6/T3	TI 1 C5/T3	TI 1 C4/T3	TI 1 C3/T3	TI 1 C2/T3	TI 1 C1/T3	PC(ANT)/T3	S/T3	Blank w/o AB
G	Blank	NC/T3	TI 1 C7/T3	TI 1 C6/T3	TI 1 C5/T3	TI 1 C4/T3	TI 1 C3/T3	TI 1 C2/T3	TI 1 C1/T3	PC(ANT)/T3	S/T3	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

523

524 Figure 5. continued

525

526 (B) the next set of plates:

527

Agonism
plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
C	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
D	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
E	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
G	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
G	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

528

Antagonism
plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	S/T3	TI 2 C7/T3	TI 2 C6/T3	TI 2 C5/T3	TI 2 C4/T3	TI 2 C3/T3	TI 2 C2/T3	TI 2 C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
C	Blank	S/T3	TI 2 C7/T3	TI 2 C6/T3	TI 2 C5/T3	TI 2 C4/T3	TI 2 C3/T3	TI 2 C2/T3	TI 2 C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
D	Blank	S/T3	TI 2 C7/T3	TI 2 C6/T3	TI 2 C5/T3	TI 2 C4/T3	TI 2 C3/T3	TI 2 C2/T3	TI 2 C1/T3	UC (PCM)	REF(T3) C1	Blank w/o AB
E	Blank	S/T3	TI 3 C7/T3	TI 3 C6/T3	TI 3 C5/T3	TI 3 C4/T3	TI 3 C3/T3	TI 3 C2/T3	TI 3 C1/T3	Ref(DPH)C1/T3	S/T3	Blank w/o AB
G	Blank	S/T3	TI 3 C7/T3	TI 3 C6/T3	TI 3 C5/T3	TI 3 C4/T3	TI 3 C3/T3	TI 3 C2/T3	TI 3 C1/T3	Ref(DPH)C1/T3	S/T3	Blank w/o AB
G	Blank	S/T3	TI 3 C7/T3	TI 3 C6/T3	TI 3 C5/T3	TI 3 C4/T3	TI 3 C3/T3	TI 3 C2/T3	TI 3 C1/T3	Ref(DPH)C1/T3	S/T3	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

529

530 **3.3.4 Measurement of cell proliferation**

531 Perform the cell proliferation assay according to Section 2: "AlamarBlue assay" in SOP
 532 "Determination of cell proliferation in T-screen assay". Before performing the assay, observe the cells
 533 under the microscope to record cytotoxic effect or precipitates (if any) as described in in SOP
 534 *Determination of cell proliferation in T-screen assay* in Section 2.1. Wells where precipitates were
 535 observed are excluded from calculation as recommended in SOP "Determination of cell
 536 *proliferation in T-screen assay*", Section 2.1.1 .

537 Observation of GH3 cell morphology after the exposure (focus on cytotoxic features, e.g. destroyed
 538 cells, destruction of the cell layers; intracytoplasmatic granules) is essential to confirm cytotoxic
 539 effect (the lower cellular activity of GH3 cells cultured in PCM medium without T3 as determined
 540 with the cell proliferation assay) and distinguish it from lack of agonistic effect (basal or low activity
 541 in PCM without T3).

542 **3.4 Data Analysis and Calculations**

543 Results of the cell proliferation assay, i.e. **% AlamarBlue reduction (%AR) or % Dye reduction (%DR)**,
 544 calculated based on optical density (OD; absorbance) for each triplicate sample using formulas given
 545 in SOP *Determination of cell proliferation in T-screen assay* (Section 2.2.1 or Section 2.2.2,
 546 respectively) are used to calculate the agonistic and antagonistic effect according to formulas given
 547 below (Section 3.4.1 and Section 3.4.2, respectively).

548 **3.4.1 Data analysis for AGONISM plate**

549 3.4.1-1 Based on received results (optical density/absorbance), calculate **% AlamarBlue reduction**
 550 **(%AR) or % Dye reduction (%DR)** for each triplicate sample using formulas given in SOP
 551 *Determination of cell proliferation in T-screen assay* in Section 2.2.1 or Section 2.2.2,
 552 respectively.

553 3.4.1-2 Calculate the increase of cell proliferation that is expressed as **the relative proliferative**
 554 **effect (RPE)** for all samples. The response observed at 2 nM T3 [Ref(T3) C1] is considered
 555 as the maximum response and set as 100%. The response for the solvent control [SC] is set
 556 at 0%. Exemplary results are presented in Figure 6A.

$$RPE = \frac{\%DR_x - \%DR_{SC}}{\%DR_{Ref(T3)C1} - \%DR_{SC}} \times 100$$

557 **OR***

$$RPE = \frac{\%AR_x - \%AR_{SC}}{\%AR_{Ref(T3)C1} - \%AR_{SC}} \times 100$$

558 where:

- 559 SC - solvent control;
 560 Ref(T3) C1 - the max. response observed for T3 (the concentration no. 1; 2 nM);
 561 x - the effect of TI, PC, NC or REF T3 at the concentration analysed

562 * - %AR or %DR is used depending on the method choose for calculation results of
 563 the proliferation assay (SOP *Determination of cell proliferation in T-screen assay* in Section
 564 2.2.1 or Section 2.2.2, respectively)

565 3.4.1-3 Determine **the EC50 value of the reference item (Ref(T3))** and **the EC50 value of the test**
 566 **item (TI)** (for dose response) or the concentration giving highest RPE value above specific
 567 limit (e.g. > 10% increase or significantly different from SC). The EC50 value can be
 568 determined using the Hill curve model in a statistic programme e.g. GraphPad. The Hill
 569 curve model is a logistic regression model (variable slope, 4 parameters) that uses the
 570 following function:

$$571 \quad y = RPE_{\min} + \frac{(RPE_{\max} - RPE_{\min})}{(1 + 10^{((\text{LogEC}_{50} - x) * \text{HillSlope}))}}$$

572 where:

- 573 x - Log of concentration
- 574 y - Relative induction (%)
- 575 RPE_{\max} - Maximum relative induction (%)
- 576 RPE_{\min} - Minimum relative induction (%)
- 577 LogEC_{50} - Log of concentration at which 50% of maximum relative induction is observed
- 578 HillSlope - Slope factor of the Hill curve

579 3.4.1-4 Calculate the Z-factor for each plate tested according the following formula:

$$580 \quad \mathbf{Z\text{-factor}}_{\text{plate no.}} = 1 - 3 \times \frac{(\text{SD}_{\text{plate no.}} [\text{SC}] + \text{SD}_{\text{plate no.}} [\text{Ref(T3)C1}])}{\text{abs} (\%AR * \text{plate no.} [\text{SC}] - \%AR * \text{plate no.} [\text{Ref(T3)C1}])}$$

582 where:

- 583 abs - absolute value;
- 584 SC - solvent control;
- 585 Ref (T3) C1 - the max. response observed for T3 (the concentration no. 1; 2 nM);
- 586 * - %AR or %DR is used depending on the method choose for calculation results of
- 587 the proliferation assay (SOP *Determination of cell proliferation in T-screen assay*
- 588 in Section 2.2.1 or Section 2.2.2, respectively)

589 3.4.2 Data analysis for ANTAGONISM plate

590 3.4.2-1 Based on received results (optical density/absorbance), calculate **% AlamarBlue reduction**
 591 **(%AR)** or **% Dye reduction (%DR)** for each triplicate sample using formulas given in SOP
 592 *Determination of cell proliferation in T-screen assay*; Section 2.2.1 or Section 2.2.2,
 593 respectively.

594 3.4.2-2 Calculate the decrease of cell proliferation expressed as **the relative inhibitory effect (RIE)**
 595 for all samples. The response for solvent control in the presence of EC50 T3 [S/T3]) is 100%;
 596 the max. response observed for REF in the ANTAGONISM plates in the presence of EC50 T3
 597 ([Ref(DPH) Cmax/T3]) is set at 0%. Exemplary results are presented in Figure 6B.

$$598 \quad \text{RIE} = \frac{\%DR_{x/T3} - \%DR_{\text{Ref(DPH)Cmax/T3}}}{\%DR_{S/T3} - \%DR_{\text{Ref(DPH)Cmax/T3}}} \times 100$$

599

600 **OR***

$$RIE = \frac{\%AR_{x/T3} - \%AR_{Ref(DPH)C_{max}/T3}}{\%AR_{S/T3} - \%AR_{Ref(DPH)C_{max}/T3}} \times 100$$

601

602 where:

603 S/T3 - solvent control in the presence of EC50 T3;

604 Ref(DPH)C_{max}/T3 - the max. response observed for the reference item (the concentration
605 no. 1; 2 nM) in the presence of EC50 T3;

606 x/T3 - the effect of TI, PC, NC or Ref(DPH) at the concentration analysed
607 in the presence of EC50 T3

608 * - %AR or %DR is used depending on the method choose for calculation
609 results of the proliferation assay (SOP *Determination of cell proliferation*
610 *in T-screen assay* in Section 2.2.1 or Section 2.2.2, respectively)

611

612 3.4.2-3 Determine the **IC50 value of the reference item (Ref(DPH))** and the **IC50 value of the test**
613 **item (TI)** (for dose response) according to the formula given below or the concentration
614 giving highest RIE value above specific limit (e.g. > 20% decrease). The IC50 value can be
615 determined using the Hill curve model in a statistic programme e.g. GraphPad. The Hill
616 curve model is a logistic regression model (variable slope, 4 parameters) that uses the
617 following function:

$$y = RIE_{min} + \frac{(RIE_{max} - RIE_{min})}{(1 + 10^{((LogIC_{50} - x) * HillSlope)})}$$

618 x = Log of concentration

619 y = Relative inhibition (%)

620 RIE_{max} = Maximum relative inhibition (%)

621 RIE_{min} = Minimum relative inhibition (%)

622 LogIC₅₀ = Log of concentration at which 50% of maximum relative inhibition is observed

623 HillSlope = Slope factor of the Hill curve

624

625 3.4.2-4 Calculate the **Z-factor** for each plate tested according the following formula.:

626

$$\mathbf{Z-factor}_{plate\ no.} = 1 - 3 \times \frac{(SD_{plate\ no.} [S/T3] + SD_{plate\ no.} [Ref(DPH) C_{max}/T3])}{abs(\%DR^*_{plate\ no.} [S/T3] - \%DR^*_{plate\ no.} [Ref(DPH)C_{max}/T3])}$$

627 where: abs - absolute value;

628 SD - Standard deviation

629 S/T3 - solvent control;

630 Ref(DPH)C_{max}/T3 - the max. response observed for the reference item (the
631 concentration no. 1; 2 nM) in the presence of EC50 T3;

632 * - %AR or %DR is used depending on the method choose for calculation results of
633 the proliferation assay (SOP *Determination of cell proliferation in T-screen assay* in
634 Section 2.2.1 or Section 2.2.2, respectively).

635
636
637

638 **3.4.3 Relative potency of test items**

639 **Agonism:**

640 *The relative potency of test items* is calculated by dividing the EC₅₀ of T3 by the EC₅₀ of the test item
 641 (or the concentration giving the highest RPE value that is above 10% increase as compared to SC).


642 **Antagonism:**

643 *The relative potency of test items* is calculated by dividing the IC₅₀ of REF by the IC₅₀ of the test item
 644 (or the concentration giving the highest RIE value that is above > 10% decrease as compared to SC).


645 **Figure 6.** Exemplary results obtained for (A) AGONISM plate 1 and (B) ANTAGONISM plate 1

646 (corrected graphs will be added after PART 1/receiving new data)

647 A



656 B



664 **3.5 Acceptance criteria**

665 **To be developed on basis of historical data with the reference and control items.**

666 **3.5.1 Acceptance criteria for AGONISM plate**

667 3.5.1-1 Mean EC₅₀ value of T3 should be -10±0.4 log₁₀(Molar) units (in the range from -10.4 to -9.6
668 log₁₀(Molar) units).

669 3.5.1-2 %DR or %AR for UC should not be more than 15% different from TI SC and REF SC

670 3.5.1-3 Z-factor > 0.5

671 **3.5.2 Acceptance criteria for ANTAGONISM plate**

672 3.5.2-1 RPE for EC₅₀ T3 should be in the range of 30-70% = will be confirmed after PART 1

673 3.5.2-2 Z-factor > 0.5

674 **4. Reference(s)**

- 675 • Ghisari, M., & Bonefeld-Jorgensen, E. C. (2005). Impact of environmental chemicals on the
676 thyroid hormone function in pituitary rat GH3 cells. *Molecular and Cellular Endocrinology*, 244:
677 31–41. <http://doi.org/https://doi.org/10.1016/j.mce.2005.01.013>
- 678 • Gutleb, A. C., Meerts, I. A. T. M., Bergsma, J. H., Schriks, M., & Murk, A. J. (2005). T-Screen as a
679 tool to identify thyroid hormone receptor active compounds. *Environmental Toxicology and*
680 *Pharmacology*, 19: 231–238. <http://doi.org/https://doi.org/10.1016/j.etap.2004.06.003>
- 681 • Sirbasku, D.A., Pakala, R., Sato, H., Eby, J.E. (1991). Thyroid hormone dependent pituitary tumor
682 cell growth in serum-free chemically defined culture. A new regulatory role for apotransferrin.
683 *Biochem.* 30: 7466–7477.
- 684 • Schriks M, Vrabie CM, Gutleb AC, Faassen EJ, Rietjens IM, Murk AJ. (2006) T-screen to quantify
685 functional potentiating, antagonistic and thyroid hormone-like activities of poly halogenated
686 aromatic hydrocarbons (PHAHs). *Toxicol In Vitro.* 20: 490-498.
- 687 • Hohenwarter O, Waltenberger A, Katinger H. (1990) An in vitro test system for thyroid hormone
688 action. *Anal. Biochem.* 234: 56-59.
- 689 • Koyanagi M, Kawakabe S, Arimura Y. (2016). A comparative study of colorimetric cell proliferation
690 assays in immune cells. *Cytotechnology*, 68, 1489-1498. doi:10.1007/s10616-015-9909-2.

GETTING IN TOUCH WITH THE EU

In person

All over the European Union there are hundreds of Europe Direct centres. You can find the address of the centre nearest you online (european-union.europa.eu/contact-eu/meet-us_en).

On the phone or in writing

Europe Direct is a service that answers your questions about the European Union. You can contact this service:

- by freephone: 00 800 6 7 8 9 10 11 (certain operators may charge for these calls),
- at the following standard number: +32 22999696,
- via the following form: european-union.europa.eu/contact-eu/write-us_en.

FINDING INFORMATION ABOUT THE EU

Online

Information about the European Union in all the official languages of the EU is available on the Europa website (european-union.europa.eu).

EU publications

You can view or order EU publications at op.europa.eu/en/publications. Multiple copies of free publications can be obtained by contacting Europe Direct or your local documentation centre (european-union.europa.eu/contact-eu/meet-us_en).

EU law and related documents

For access to legal information from the EU, including all EU law since 1951 in all the official language versions, go to EUR-Lex (eur-lex.europa.eu).

Open data from the EU

The portal data.europa.eu provides access to open datasets from the EU institutions, bodies and agencies. These can be downloaded and reused for free, for both commercial and non-commercial purposes. The portal also provides access to a wealth of datasets from European countries.

The European Commission's science and knowledge service

Joint Research Centre

JRC Mission

As the science and knowledge service of the European Commission, the Joint Research Centre's mission is to support EU policies with independent evidence throughout the whole policy cycle.



EU Science Hub
joint-research-centre.ec.europa.eu



@EU_ScienceHub



EU Science Hub - Joint Research Centre



EU Science, Research and Innovation



EU Science Hub



EU Science