


STANDARD OPERATING PROCEDURE

T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin, version 2.0

EURL ECVAM validation study of a battery of mechanistic methods relevant for the detection of chemicals that can disrupt the thyroid hormone system

Roszak J.



This Standard Operating Procedure (SOP) has been prepared within the context of a collaboration agreement with the Joint Research Centre (JRC) Directorate for Health, Consumers and Reference Materials (Chemicals Safety and Alternative Methods Unit F3 / EURL ECVAM), for the validation of mechanistic methods to identify potential modulators of thyroid hormone signalling. It aims to provide evidence-based scientific support to the European policymaking process. The contents of this publication do not necessarily reflect the position or opinion of the European Commission. Neither the European Commission nor any person acting on behalf of the Commission is responsible for the use that might be made of this publication. For information on the methodology and quality underlying the data used in this publication for which the source is neither Eurostat nor other Commission services, users should contact the referenced source. The designations employed and the presentation of material on the maps do not imply the expression of any opinion whatsoever on the part of the European Union concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries.

This SOP is part of a series of 3 SOPs used to perform the "T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin":

1. **SOP "T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin" v2.0** (used in Part 2 of the validation study)
2. SOP "Handling, Maintenance and Quality Control of the GH3 cell line " v2.0 (used in Part 2 of the validation study)
3. SOP "Determination of cell proliferation in T screen assay" v2.0 (used in Part 2 of the validation study)

The method was developed by Arno Gutleb, LIST (Luxembourg) and subsequently implemented by the EU-NETVAL test facility NIOM (Poland) within the validation study. For the purpose of running the method on a robotic platform during part 2 of the validation study, EURL ECVAM prepared a 4th SOP.

4. SOP "T-Screen assay using a QUANTITATIVE HIGH THROUGHPUT screening approach" v1.0 (used in Part 2 of the validation study)

Contact information

Method developer:

Arno Gutleb
PhD, ERT
Luxembourg Institute of Science and Technology (LIST)
Environmental Research and Innovation (ERIN) Department
41, rue du Brill
L-4422 Belvaux
Grand-duchy of Luxembourg
Tel: +352 275 888 481
e-mail: arno.gutleb@list.lu

EU-NETVAL facility

Joanna Roszak
National Center for Alternative Methods to Toxicity Assessment (CMA)
Department of Translational Research
Nofer Institute of Occupational Medicine
91-348 Lodz, 8 Sw. Teresy St, POLAND
e-mail: Joanna.Rozsak@imp.lodz.pl
<https://www.imp.lodz.pl/krajowe-centrum-metod-alternatywnych-do-oceny-toksycznosci>

EU Science Hub

<https://joint-research-centre.ec.europa.eu>

JRC133177

Ispra: European Commission, 2023
© European Union, 2023



The reuse policy of the European Commission documents is implemented by the Commission Decision 2011/833/EU of 12 December 2011 on the reuse of Commission documents (OJ L 330, 14.12.2011, p. 39). Unless otherwise noted, the reuse of this document is authorised under the Creative Commons Attribution 4.0 International (CC BY 4.0) licence (<https://creativecommons.org/licenses/by/4.0/>). This means that reuse is allowed provided appropriate credit is given and any changes are indicated.

For any use or reproduction of photos or other material that is not owned by the European Union, permission must be sought directly from the copyright holders. The European Union does not own the copyright in relation to the following elements:
- Cover page illustration, © BioRender.com

How to cite this report: Roszak J., Standard Operating Procedure for T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin, version 2.0, applied in Part 2 of the EURL ECVAM thyroid validation study, European Commission, Ispra, 2023, JRC133177.

1
2
3
4
5
6
7
8
9
10
11
12
13
14


Standard Operating Procedure

T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin

Author: EU-NETVAL laboratory NIOM, Joanna Roszak

Version: 02

Date: 25.01.2022



CONFIDENTIAL

15 **Contents**

16 1. Description of the method 4

17 1.1 Scope and Limitations of the Application 4

18 1.2 Safety 4

19 1.3 Required Standard Operating Procedures 4

20 1.4 Test System 4

21 1.5 Apparatus and Materials 4

22 1.6 Reagents and Chemicals 5

23 1.7 Reference items, positive and negative control Item(s) 5

24 1.8 Preparations of media and reagents 5

25 1.8.1 Cell culture medium (cDMEM/F12; completed cDMEM/F12) 5

26 1.8.2 PCM medium 6

27 1.8.3 Stock and working solutions of chemicals 6

28 1.8.4 Stock solutions of reagents 7

29 2. Verification of the EC50 value of T3 8

30 2.1 Seeding the GH3 cells onto 96-well plate 8

31 2.2 Exposure of the cells to the range of T3 concentrations 8

32 2.3 Measurement of cell proliferation 8

33 2.4 Acceptance criteria 9

34 3. T-Screen experimental procedure 9

35 3.1 Preparation of Test Item 9

36 3.1.1 Stock solutions 9

37 3.1.2 Work solutions to determine solubility 10

38 3.2 Pre-screen experiment 10

39 3.2.1 Seeding the GH3 cells onto 96-well plate 10

40 3.2.2 Preparation of test, reference and control item solutions 11

41 3.2.3 Exposure of the cells 11

42 3.2.4 Measurement of cell proliferation and data analysis 12

43 3.2.5 Acceptance criteria for pre-screen assay 12

44 3.2.6 Identification of the range concentrations of TI for the next experiments 12

45 3.3 T-Screen assay 15

46 3.3.1 Seeding the GH3 cells onto 96-well plate 15

47 3.3.2 Preparation of test, reference and control item solutions 15

48 3.3.3 Exposure of the cells 15

49 3.3.4 Measurement of cell proliferation 19

50 3.4 Data Analysis and Calculations 19

51 3.5 Acceptance criteria 21

52 4. Reference(s) 21

53 5. A change log 22

54

55 *Figure 1. Scheme of the plate layout for verification of the EC50 value of T3 7*

56 *Figure 2. The T-screen assay workflow with absorbance of AlamarBlue dye detection. 8*

57 *Figure 3. A schema of test item preparation for pre-screen experiment and the T-Screen test. 12*

58 *Figure 4. Plate layout for pre-screen experiment 13*

59 *Figure 5. Plate layout for the 8a T-Screen assay 16*

60 *Figure 6. Exemplary results obtained for reference item with agonistic potential 19*

61 **1. Description of the method**

62 The T-Screen represents an *in vitro* bioassay based on thyroid hormone (TH) dependent cell
63 proliferation of a rat pituitary tumour cell line (GH3) in serum-free medium. It can be used to study
64 interference of compounds with TH at the cellular level, thus bridging the gap between limitations of
65 assays using either isolated molecules (enzymes, transport proteins) or complex *in vivo* experiments
66 with all the complex feedback mechanisms present. Test items are tested in the absence of TH to
67 test for agonistic potency.

68 GH3 cell growth is increased in the presence of TH agonists. Cell growth is measured with
69 AlamarBlue/Resazurin cell proliferation assay using a standard plate reader. In this method a
70 colorimetric assay is used, where resazurin is reduced from a blue oxidized form into its violet
71 reduced form of resorufin. The change of colour can be detected as a change in absorbance using a
72 microplate reader.

73 **1.1 Scope and Limitations of the Application**

74 The T-screen assay is used for testing of compounds for TH receptor (TR)-mediated effects.

75 **1.2 Safety**

76 The use of endocrine disrupting chemicals can be extremely hazardous, and precautions such as
77 using gloves, protective goggles and masks under a laminar flow hood should always be taken while
78 performing chemical treatments.

79 **1.3 Required Standard Operating Procedures**

- 80 • *Handling and Maintenance of GH3 cell line*
- 81 • *Determination of cell proliferation in T-screen assay*

82 **1.4 Test System**

- 83 • GH3-cell line – a rat pituitary tumour cell line, (ATCC® CCL-82.1™)

84

85 The GH3 cells may be used when the following requirements are met:

86 **1.4-1** Test system must be free of microbial and mycoplasma contamination

87 **1.4-2** The doubling time of the GH3 cells in cDMEM/F12 medium should be 42 ± 5 h

88 **1.4-3** Cell number in PCM after 72 ± 1 h and 96 ± 1 h of culture should be at least 35% lower
89 that cell number in cDMEM/F12, when determined simultaneously.

90 **1.4-4** Cells should be used in passage from 3 to 25

91 **1.5 Apparatus and Materials**

- 92 • Balance (minimal value: 50 mg; the precision requirement: $50 \pm \leq 0.001$ mg) and the
93 appropriate mass standards F1, e.g. 1 - 500 mg $\pm \leq 0.02$ mg; 1 g $\pm \leq 0.01$ mg
- 94 • Clean glass vials for the preparation of stock solutions (e.g., gas chromatography vials or
95 Sarstedt #86.1509)
- 96 • 96-well microplates (e.g., Nunc #167008)
- 97 • Freezer below -16 °C
- 98 • Refrigerator at $2-10$ °C
- 99 • CO₂ humidified incubator at $37^{\circ}\text{C} \pm 2$ °C, 5% CO₂ $\pm 0.5\%$
- 100 • 37 °C water bath
- 101 • PipetteAid
- 102 • Pipettors (p1000, p200, p20) or Micropipettes (1000-100; 100-10; 10-0.5)
- 103 • Centrifuge
- 104 • Vacuum aspirator

- 105 • Microplate Reader; for measuring absorbance
- 106 • Laminar Flow Hood
- 107 • Culture Flasks (T75, T25; e.g. Nunc, Falcon)
- 108 • Serological pipettes (1 mL, 5 mL, 10 mL, and 25 mL)
- 109 • Sterile, filter pipette tips (10-20 µL, 200 µL, and 1000 µL)
- 110 • Sterile, pipette tips without filter (e.g. 10-20 µL, 200 µL, and 1000 µL)
- 111 • Conical tubes (e.g. 5 mL, 15 mL, and 50 mL)
- 112 • Polypropylene Cluster Tubes (1.1 mL e.g. Corning #MTS-11-8-C)
- 113 • Syringe filters (0.22µm)

114 1.6 Reagents and Chemicals

- 115 • Foetal calf/bovine serum (e.g. Gibco #10270106)
- 116 • Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12; 1:1 mixture
- 117 with phenol red and HEPES (15 mM) (e.g. Sigma-Aldrich #D6421)
- 118 • DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES and w/o phenol red (e.g. Sigma-Aldrich
- 119 #D6434)
- 120 • HEPES 1M solution (e.g. Sigma-Aldrich #H0887)
- 121 • Penicillin-Streptomycin Solution (e.g. Sigma-Aldrich #P0781)
- 122 • Bovine insulin (e.g. Sigma-Aldrich #I0516)
- 123 • Ethanolamine (e.g. Sigma-Aldrich #E0135)
- 124 • Sodium selenite (e.g. Sigma-Aldrich #S5261)
- 125 • Human apotransferrin (e.g. Sigma-Aldrich #T2036)
- 126 • Bovine serum albumin (e.g. Sigma-Aldrich #A9418)
- 127 • DMSO 99.9% purity (e.g. Sigma-Aldrich #D8418)
- 128 • Ethyl alcohol 99.8% pure p.a. (e.g. POCH #396480111)
- 129 • Deionised water (dH₂O)

130 1.7 Reference items, positive and negative control Item(s)

	Agonism
Reference item	3,3'-5-triiodothyronine (T3) [CAS 6893-02-3], ≥ 95% purity
Positive control item	T4, L-Thyroxine (3,3',5,5''-Tetraiodo-L-thyronine) [CAS 51-48-9], ≥ 98%
Negative control item	Furosemide (Fs) [54-31-9]
Cytotoxic positive control item	Sodium Dodecyl Sulfate (SDS) [CAS 151-21-3]

132

133 1.8 Preparations of media and reagents

134 1.8.1 Cell culture medium (cDMEM/F12; completed cDMEM/F12)

- 135 • Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12; 1:1
- 136 mixture; with phenol red and 15 mM HEPES) supplemented with 10% Foetal calf serum, 2.5
- 137 mM L-Glutamine and additional amount of HEPES to obtain 25 mM.

138

139 To prepare *cDMEM/F-12* the following volume of supplements should be added:

	Final concentration	Volume [mL]
DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES (e.g. Sigma-Aldrich #D6421)*		500
Heat inactivated FCS/FBS**	10%	57.6
200 mM L-Glutamine	2.5 mM	7.2
1 M HEPES	up to 25 mM	5.8
Penicillin-Streptomycin solution (100x)	1%	5.8
Total		576.4

140 *if different DMEM/F-12 is used the final concentration of HEPES and L-Glutamine should be
141 adjusted to 2.5 mM and 25 mM, respectively.

142 ** if heat-inactivated FCS/FBS was purchased, it should be defrost at 2-10°C e.g. during the night and
143 stored in aliquots at below -16°C.

144 If non-inactivated FCS/FBS was purchased, it should be heat-inactivated in water bath at 56-57°C for
145 30-35 min and filtered (0.2 µm) before aliquoting.

146

147 1.8.2 PCM medium

148 PCM medium is a Dulbecco's modified Eagle's medium/Ham's Nutrient Mixture F-12 (DMEM/F-12;
149 1:1 mixture) without phenol red supplemented with 10 µg/mL bovine insulin, 10 µM ethanolamine,
150 10 ng/mL sodium selenite, 10 µg/mL human apotransferrin, 500 µg/mL bovine serum albumin.
151 Neither T3 nor T4 are added to the PCM medium.

152 To prepare *PCM medium* the following volume of supplements should be added:

	Final concentration	Volume [mL]
DMEM/F-12 w/o L-Glutamine, with 15 mM HEPES (e.g. Sigma-Aldrich # D6434)*		500
200 mM L-Glutamine	2.5 mM	6.5
1 M HEPES	up to 25 mM	5.2
Penicillin-Streptomycin solution (100x)	1%	5.2
10 mg/mL bovine insulin	10 µg/ml	0.523
50 mM ethanolamine (Section 1.8.4-1)	10 µM	0.105
50 µg/mL sodium selenite (Section 1.8.4-2)	10 ng/mL	0.105
2 mg/mL human apotransferrin (Section 1.8.4-3)	10 µg/mL	2.6
100 mg/mL bovine serum albumin (Section 1.8.4-4)	500 µg/mL	2.6
Total		523.0

153 *if different DMEM/F-12 is used the final concentration of HEPES and L-Glutamine should be
154 adjusted to 2.5 mM and 25 mM, respectively.

155 1.8.3 Stock and working solutions of chemicals

156 To prepare the **exposure concentration** of every chemical given below, firstly the **concentrated**
157 **solution** is prepared that is used to prepare **the working solution** (the double desired exposure
158 concentration) according to Section 3.1.1 and/or Section 3.1.2.

- 159 1.8.3-1 **3,3'-5-triiodothyronine (T3) [2mM stock solution] and [4 nM working solution]**
160 Dissolve T3 in DMSO to produce a 1 mM stock solution. Use fresh or store in aliquots below
161 -16 °C.
162 T3 is used in the T-Screen assay in the full dose range (the first Agonism plate) and in the
163 highest test concentration, i.e. 2 nM – to prepare the concentrated solution of T3 [1 µM],
164 dilute T3 stock solution in the solvent/DMSO as follows: 2 mM → (100x) → 20 µM --
165 (20x)→ 1 µM. Then, use the *concentrated solution* of T3 to prepare *working solution*
166 according to Section 3.2.2
- 167 1.8.3-2 **3,3',5,5'-tetraiodothyroxine (T4) [2 mM stock solution] and [20 nM working solution]**
168 Dissolve T4 in DMSO to produce a 2 mM stock solution. Use fresh or store in aliquots below
169 -16 °C.
170 The exposure concentration of T4 in the T-Screen assay is 10 nM – to prepare the
171 concentrated solution of T4 [5 µM], dilute T4 stock solution in the solvent/DMSO as
172 follows: 2 mM → (400x) → 5 µM. Then, use *the concentrated solution* of T4 to prepare
173 *working solution* according to Section 3.2.2.
- 174 1.8.3-3 **Furosemide (FS) [100 mM stock solution] and [200 nM working solution]**
175 Dissolve FS in DMSO to produce a 100 mM stock solution. Use fresh or store in aliquots
176 below -16 °C.
177 The exposure concentration of FS in the T-Screen assay is 100 nM – to prepare the
178 concentrated solution of Fs [50 µM], dilute Fs stock solution in the solvent/DMSO as
179 follows: 100 mM → (100x)→ 1 mM → (20x)→ 50 µM. Then, use *the concentrated solution of*
180 *Fs* to prepare *working solution* according to Section 3.2.2.
- 181 1.8.3-4 **Sodium Dodecyl Sulfate (SDS) [100 mM stock solution] and [200 µM working solution]**
182 Prepare the 100 mM stock solution of SDS in DMSO. Use fresh or store in aliquots at room
183 temperature.
184 The exposure concentration of SDS in the T-Screen assay is 100 µM – to prepare the
185 concentrated solution of SDS [50 mM], dilute SDS stock solution in the solvent/DMSO as
186 follows: 100 mM → (2x)→ 50 mM. Then, use the *concentrated solution* of SDS to prepare
187 *working solution* according to Section 3.2.2
- 188 **1.8.4 Stock solutions of reagents**
- 189 1.8.4-1 **Ethanolamine [50 mM]**
190 Prepare a 50 mM stock solution of ethanolamine in dH₂O by diluting a pure ethanolamine
191 (16.6 M) 332x, i.e. 10 µL of ethanolamine (16.6 M) added to 3.310 mL of dH₂O, filter (0.22
192 µm pore size). Use fresh or store at room temperature in closed glass containers.
- 193 1.8.4-2 **Sodium selenite [50 µg/mL]**
194 Prepare a stock solution of sodium selenite (50 µg/mL) in dH₂O, filter (0.22 µm pore size)
195 and store in aliquots below -16 °C.
- 196 1.8.4-3 **Human apotransferrin [2 mg/mL]**
197 Prepare a 2 mg/mL stock solution of human apotransferrin in dH₂O, filter (0.22 µm pore
198 size), and store in aliquots below -16 °C. Stock solutions of apo-transferrin stored at 2-8°C
199 are stable for 5–10 days.
- 200 1.8.4-4 **Bovine serum albumin [100 mg/mL]**
201 Prepare a 100 mg/mL stock solution of bovine serum albumin in dH₂O, filtered (0.22 µm
202 pore size) and store in aliquots at 2-8°C

203 **2. Verification of the EC50 value of T3**

204 Once for each working cell bank and for each new batch of T3 the laboratory should verify the EC50
205 value for reference item T3 as described below.

206 **2.1 Seeding the GH3 cells onto 96-well plate**

207 Perform the cell seeding as described for the T-screen (Section 3.3.1). Seed cells onto a 96-well
208 microplate in three replicates (e.g. rows B-D; Figure 1).

209 **2.2 Exposure of the cells to the range of T3 concentrations**

210 2.2 - 1. Prepare T3 at the range of double desired concentrations by diluting the concentration
211 of 4 nM using dilution factor (DF) = 3. The final 7 concentration of T3 should be: 2;
212 0.667; 0.222; 0.074; 0.025; 0.008 and 0.003 nM

213 2.2 - 2. Add 100 µL T3 in PCM medium at the appropriate concentration into the 100 µL already
214 present (with the cells) in the respective well(s) (to achieve a 1:1 dilution). Test T3 and
215 control items in triplicate (Figure 1).

216 2.2 - 3. Incubate the plate for 96 ± 2 h at 37 ± 1 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid
217 atmosphere.
218

219 **Figure 1. Scheme of the plate layout for verification of the EC50 value of T3**

	1	2	3	4	5	6	7	8	9	10	11	12
A	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z
B	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
C	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
D	Blank	SC	T3 C7	T3 C6	T3 C5	T3 C4	T3 C3	T3 C2	T3 C1	PCM	SC	Blank w/o AB
E	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z	Z
F												
G												
H												

220 **SC** (wells 2B-D and 11B-D) - solvent control

221 **PCM** (wells 10B-D) - test system control

222 **T3 C** (wells 3B-9D) - the range of T3 concentrations (0.003 - 2 nM); where C1 is the lowest T3
223 concentration

224 **Z** – external wells (medium or PBS)

225 **Blank** - PCM medium without cells (AlamarBlue is added in the proliferation assay); used for
226 calculations of % Dye reduction (%DR) according to formulas given in SOP *Determination of cell*
227 *proliferation in T-screen assay* in Section 2.

228 **Blank w/o AB** - PCM medium without cells (AlamarBlue is not added in the proliferation assay);
229 used for calculations of %DR according to formula given in SOP *Determination of cell*
230 *proliferation in T-screen assay* in Section 2.

231 **2.3 Measurement of cell proliferation**

232 2.3 - 1. Perform the cell proliferation assay according to Section 2: “AlamarBlue assay”
233 described in SOP *Determination of cell proliferation in T-screen assay*”.

234 2.3 - 2. Based on received results (optical density/absorbance), calculate **% Dye reduction**
 235 **(%DR)** for each triplicate sample using formulas given in SOP *Determination of cell*
 236 *proliferation in T-screen assay* in Section 2.

237 2.3 - 3. Calculate the increase of *Cell proliferation* that is expressed as **the relative proliferative**
 238 **effect (RPE)** according to the formula given in Section 3.4.2

239 2.3 - 4. Calculate the **EC50** of T3 according to the formula given in Section 3.4.3.

240 2.4 Acceptance criteria

241 Based on the obtained results the mean EC50 value of T3 should be calculated. The mean EC50
 242 value should be $-10 \pm 0.4 \log_{10}(\text{Molar})$ units (in the range from -10.4 to -9.6 $\log_{10}(\text{Molar})$ units).

243 If acceptance criteria for the mean EC50 value of T3 are met, 0.1 nM T3 should be used as the
 244 EC50 value of T3 in the T-screen test.


245 Every result, i.e. the EC50 of T3, should be registered to create the historical data EC50 value of T3.

246 Also, the mean EC50 of T3 calculated from the results obtained during every T-screen test (the first
 247 plate, Section 3.4.3) should be included to the historical data.

248 3. T-Screen experimental procedure

249 Figure 2 summarises the main steps required to run this method.

250 **Figure 2.** The T-screen assay workflow with absorbance of AlamarBlue dye detection.



251

252

253 3.1 Preparation of Test Item

254 The test item's solubility in the solvent and medium will need to be assessed. Both stock solutions
 255 and working solutions should be assessed.

256 3.1.1 Stock solutions

257 *It is recommended to use the same solvent for all reference and control items and to use that solvent*
 258 *for the test items, to ensure all samples are tested under the same conditions. In case another solvent*
 259 *is used for the test item, the impact on the test system and the results must be assessed.*

260 3.1.1-1 Prepare a 100 mM or 10 mg/mL stock solution (may be adapted on basis of data if too low
 261 or high) of test item by weighing a nominal amount into a clear glass vial and add the
 262 appropriate volume of solvent (use DMSO as a default solvent and if not possible to

263 dissolve test item in DMSO then other solvents, e.g. water, PCM medium or ethanol
264 should be used. For visual inspection, it is important to have at least 0.5 ml in the vial to be
265 able to observe well.

266 3.1.1-2 Visually check if dissolved. If not dissolved, heat (37°C; up to 30 min) may be applied to aid
267 solubility.

268 3.1.1-3 If not possible to solubilise, prepare a lower (e.g. 1:2) stock solution (or dilute existing
269 stock) and again check solubility. Change solvent if needed.

270 3.1.1-4 Continue until the stock solution is soluble.

271 **3.1.2 Work solutions to determine solubility**

272 3.1.2-1 Dilute the stock solution of test item in solvent prepared in Section 3.1.1. Prepare the
273 range of 4 concentrations (*the concentrated solutions*) using a dilution factor 2 (50 µL of
274 stock solution of the appropriate *concentrated solution* to 50 µL solvent in a 96-well plate).
275 For test item dissolved in DMSO or ethanol prepare the 500x *concentrated solutions*,
276 whereas for test item dissolved in water or PBS, prepare 100x *concentrated solutions*.

277 3.1.2-2 Prepare the 2x *concentrated solutions (working solutions)* in PCM medium by transferring:
278 2 µL of the appropriate *concentrated solutions* to 498 µL of PCM in a 24-well plate or 1.1
279 mL cluster tubes.

280 Solvent concentration in the *working solutions* of test item should not exceed 0.4% for
281 DMSO and ethanol or 10% for water or PBS (the final concentration 0.2% or 5%,
282 respectively for DMSO/ethanol or water/PBS).

283 3.1.2-3 Visually check (using a microscope) if the working solutions are dissolved.

284 3.1.2-4 Incubate the plate with the working solutions for 24 h at 37 ± 1 °C and $5 \pm 0.5\%$ (v/v) CO₂
285 in a humid atmosphere to check if any precipitates appear.

286 3.1.2-5 Identify the highest soluble concentration.

287 **3.2 Pre-screen experiment**

288 During the pre-screen experiment it is assessed if the test item is:

- 289 • A Thyroid Hormone Agonist
- 290 • Cytotoxic

291 and at which dose-range it increases or decreases cell proliferation (Range finder).
292

293 **Note:** GH3 cells are not able to divide properly in PCM without T3 but keep basal or low
294 activity. Concentrations of test items lower than 15% of the cellular activity of GH3 cells
295 cultured in PCM medium without T3 (as determined with the cell proliferation assay) are
296 considered to be cytotoxic.
297

298 **3.2.1 Seeding the GH3 cells onto 96-well plate**

299 Perform the cell seeding as described for the T-screen (Section 3.3.1).

- 300 - Seed cells into all internal wells, i.e. B2-G11.
- 301 - Add 100 µL of PCM medium to each well in columns 1 and 12, for *Blank* and *Blank w/o AB*
302 (*Figure 4*) that are used for calculations in the proliferation assay.
- 303 - Add 100 µL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
304 proper humidity for the cells.

305 **3.2.2 Preparation of test, reference and control item solutions**

306 On the day of treatment, prepare the range of 7 *working solutions*:

307 3.2.2-1 Stage 1: Prepare stock solution of test item at the concentration and solvent determined
308 above (Section 3.1.1). Prepare a dilution series of test item stock in the solvent (the range
309 of 7 *concentrated solutions*) using a dilution factor 10 (10 µL of stock solution of the
310 appropriate concentrated solution to 90 µL solvent in a 96-well plate).

311 For test item dissolved in DMSO or ethanol prepare the 500x concentrated solutions,
312 whereas for test item dissolved in water or PBS, prepare 100x concentrated solutions.

313 In the first experiment a dose range finding is performed to capture the whole dose response range
314 for the next experiments. The highest concentration to be tested is the highest non cytotoxic soluble
315 concentration and dilution factor 10 is applied for the endpoint measurement. In the next
316 experiment a dose range includes 7 concentrations where the highest concentration shows the max
317 effect (the induction of proliferation) and the lowest concentration shows no effect (effect
318 comparable to solvent control). Dilution factor is adjusted to cover 7 concentrations and each test
319 item have its own dilution factor (e.g. DF 1.5, 2, 3, 4 or 5).

320
321 3.2.2-2 Stage 2: Prepare *working solutions* as described in Section 3.1.2

322 In the first experiment/the dose range finding an interference of test item with the
323 assay/AlamarBlue reagent should be assessed. To this end, prepare *working solutions* as described in
324 Section 3.1.2 using double volumes to have enough solutions both for test plate and the additional
325 plate for testing interference of AlamarBlue with test item.

- 326
327 - Add 498 µL PCM medium into 1.1 mL-cluster tubes, then
328 - Add 2 µL of the appropriate *concentrated solution* to the mixture of PCM medium and PCM-
329 T3 medium already present into 1.1 mL-cluster tubes. Solvent concentration should not
330 exceed 0.5%.

331

332 **3.2.3 Exposure of the cells**

333 Test all samples in triplicate. Assess one test item or reference item on an upper part of a
334 plate(B2:D11) and another test item on a lower part of a plate (E2:G11). Perform the exposure as
335 follows:

336 3.2.3-1 Test plate:

- 337 - Add 100 µL *working solutions* (the double desired exposure concentrations) of test item
338 (TI) or solvent control (SC) prepared in PCM medium to the 100 µL already present (with
339 the cells) in the respective well(s) (for TI 1/SC 3B:9D/2B-D and 11B-D and for and TI 2/SC
340 3E:9G/2E:G; to achieve a 1:1 dilution).
341 - Add 100 µL PCM medium only to three wells designed "UC (PCM)" (10E:G; untreated
342 control/cell system control).
343 - Add 100 µL 2 x EC50T3 = the concentration 0.2 nM T3 (if it was positively verified; Section
344 2.4)
345 - Add 100 µL 2x Ref(T3)C1 = the highest concentration in the range described in Section
346 3.2.2 (Figure 4).

347 - Add 100 µL of PCM medium to the 100 µL already present (without the cells) in each well
348 in columns 1 and 12 (Blank and Blank w/o AB; Figure 4) that are used for calculations in the
349 proliferation assay.

350 -Add 100 µL of PCM medium or PBS to each well in rows A and H (external wells) to ensure
351 proper humidity for the cells.

352 3.2.3-2 Additional plate (interference assessment)

353 Prepare the additional plate for testing interference of AlamarBlue reagent with test items
354 (samples incubated in PCM medium without cells) – add 100 µL of PCM medium into every
355 well needed (triplicates/sample). Then, add all prepared samples (100 µL) in the same
356 order as indicated in Section 3.2.3-1. Also, add extra PCM medium in columns 1 and 12 and
357 in rows A and H, as given in Section 3.2.3-1.

358 3.2.3-3 Incubate plates for 96 ± 2 h at 37 ± 1 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid atmosphere.

359 *Plate layout for Pre-screen experiment is presented in Figure 4. It is advised to assess each item in*
360 *triplicates.*

361 **3.2.4 Measurement of cell proliferation and data analysis**

362 3.2.4-1 Perform the AlamarBlue assay according to Section 2 in SOP "Determination of cell
363 proliferation in T-screen assay" to assess cell proliferation or cytotoxic effects.

364 3.2.4-2 Based on received results (optical density/absorbance), calculate % **Dye reduction (%DR)** for
365 each triplicate sample using formulas given in SOP *Determination of cell proliferation in T-*
366 *screen assay* in Section 2.

367 3.2.4-4 Calculate the increase of cell proliferation that is expressed as **the relative proliferative**
368 **effect (RPE)** for "EC50 T3" according to Section 3.4.2

370 **3.2.5 Acceptance criteria for pre-screen assay**

371 3.2.5-1 %DR for TI SC should not be more than 15% different from UC (PCM)

372 3.2.5-2 RPE for EC50 T3 (i.e. 0.1 nM T3) should be in the range of 30-70%

374 **3.2.6 Identification of the range concentrations of TI for the next experiments**

375 Identify the range concentrations of TI for the next experiments that captures the whole dose
376 response range (from no effect to the highest agonistic effect). The range of concentration should
377 have not more than two concentrations given the maximal response and not more than two
378 concentrations given the minimal response (near background/ not more than 5% of SC).

379 3.2.6-1 If a strong cytotoxic effect is detected, i.e. three or more concentrations are found to be
380 cytotoxic, repeat the test using the changed range of concentrations (start from the highest
381 non-cytotoxic concentration and adjust DF to get 7 concentrations) to better identify range
382 of test concentrations of test item. The highest concentration for the T-screen assay should
383 be the highest noncytotoxic concentration of test item.

384 3.2.6-2 If no cytotoxic effect is observed, but the agonistic effect is not observed choose for the next
385 experiment the same range of concentration and DF 10


386 3.2.6-3 If no cytotoxic effect is observed, but any agonistic effect is observed, choose DF to obtain
387 the whole dose response (e.g. DF 1.5, 2, 3, 4 or 5) (if any).

388

389 **Figure 3.** A schema of test item preparation for pre-screen experiment and the T-Screen test.


390

Stage 1. Dilution of a stock solution to prepare the range of **the concentrated solutions**



The solvent concentration in the diluted stock solutions is considered as 100%.

Stage 2. Preparation of the **working solutions** (the double desired exposure concentrations) Cells should be exposed to a test item from minimal to maximal concentration.



The solvent concentration in **the working solutions** is considered as 0.4 %.

391 **Figure 4.** Plate layout for pre-screen experiment

	1	2	3	4	5	6	7	8	9	10	11	12	
A	*	*	*	*	*	*	*	*	*	*	*	*	
TI 1	B	Blank	SC	TI1 C7	TI1 C6	TI1 C5	TI1 C4	TI1 C3	TI1 C2	TI1 C1	UC (PCM)	SC	Blank w/o AB
	C	Blank	SC	TI1 C7	TI1 C6	TI1 C5	TI1 C4	TI1 C3	TI1 C2	TI1 C1	UC (PCM)	SC	Blank w/o AB
	D	Blank	SC	TI1 C7	TI1 C6	TI1 C5	TI1 C4	TI1 C3	TI1 C2	TI1 C1	UC (PCM)	SC	Blank w/o AB
TI 2	E	Blank	SC	TI2 C7	TI2 C6	TI2 C5	TI2 C4	TI2 C3	TI2 C2	TI2 C1	EC50 T3	Ref(T3)C1	Blank w/o AB
	F	Blank	SC	TI2 C7	TI2 C6	TI2 C5	TI2 C4	TI2 C3	TI2 C2	TI2 C1	EC50 T3	Ref(T3)C1	Blank w/o AB
	G	Blank	SC	TI2 C7	TI2 C6	TI2 C5	TI2 C4	TI2 C3	TI2 C2	TI2 C1	EC50 T3	Ref(T3)C1	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*	

392

UC (PCM)	Untreated Control; Tested as triplicate sample on each plate, PCM Medium Only
EC50 T3	PCM medium + EC50 T3; Tested as triplicate sample on each plate
SC	PCM medium + solvent; Tested as double triplicate sample on each side of the plate
Ref(T3)C1	the max. concentration of T3 used in the study/the concentration no. 1; 2 nM); it will be used for calculation of RPE for EC50
*	Outer wells with PBS or PCM medium only; Outer wells may only be used only when the plate has the additional collar/space to fill it with liquid to ensure the cells proper humidity
TI1 C [7-1]	Test item in the range of concentration (C) from C7/Cmin to C1/Cmax (C7 is the lowest concentration tested)
TI2 C [7-1]	Another test item in the range of concentration (C) from C7/Cmin to C1/Cmax (C7 is the lowest concentration tested)
Blank	PCM medium without cells (AlamarBlue is added in the proliferation assay); used for calculations of % DR according to formulas given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.
Blank w/o AB	PCM medium without cells (AlamarBlue is not added in the proliferation assay); used for calculations of %DR according to formula given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.

393 **3.3 T-Screen assay**

394 IMPORTANT! Complete Reference item standard curve (Ref(T3)) in the should be included in each
395 assay.

396 **3.3.1 Seeding the GH3 cells onto 96-well plate**

397 Number of plates depends on number of test items. Because complete Reference item standard
398 curve should be included in each assay thus the first plate is used to assess one TI together with REF
399 (T3). Another plate is used to assess the next two TI (Figure 5).

400 3.3.1-1 48 ± 2 hours prior to plating the cells onto 96-well microplates for the experiment, change
401 the standard culture medium to PCM medium (Section 1.8.2).

402 3.3.1-2 Release cells as described in SOP "*Handling and Maintenance of GH3 cell line*", Section
403 2.3.4 "Detachment of GH3 cells". Pipet cell suspension very carefully several times (cells
404 easily detach but are sensitive to shaking as well as the effects of trypsin). During
405 detachment use only PCM medium (do not use culture medium with serum).

406 3.3.1-3 Determine density of cell suspension taking into account viable cells as described in SOP
407 "*Handling and Maintenance of GH3 cell line*", Section 2.3.6.

408 3.3.1-4 Calculate cell viability as described in SOP "*Handling and Maintenance of GH3 cell line*",
409 Section 2.3.6.

410 3.3.1-5 Only if viability of GH3 is more than 90%, dilute cell suspension in PCM medium and seed
411 cells onto a 96-well microplate at a density of 2500 viable cells/well in 100 µL PCM
412 medium.

413 3.3.1-6 Pre-incubate for 24 hours ± 2 hours at 37 ± 1 °C and 5 ± 0.5% (v/v) CO₂ in a humid
414 atmosphere to allow cells to attach to bottom of wells before the treatment.

415 3.3.1-7 Before the treatment, based on microscopic observation, the following acceptance
416 criterion should be met (otherwise, the plate is rejected): not more than 50% of cells in the
417 well are floating

418 **3.3.2 Preparation of test, reference and control item solutions**

419 Prepare test item *working solutions* as for the pre-screen experiment, using the dilution factor (DF)
420 identified in the range finder pre-screen experiment that will capture the whole dose response. If no
421 agonistic effect is observed the same range of concentrations (and the same DF) should be used in T-
422 screen assay. If any agonistic effect is observed adjust DF to obtain the whole dose response (e.g. DF
423 2, 3, 4 or 5) (if any) and then repeat the experiment.

424 Reference and control items are used in T-Screen assay using the fixed concentration – the
425 preparation of *working solutions* of reference and control items (*the double desired exposure*
426 *concentrations*) are described in Section 1.8.3.

427 **3.3.3 Exposure of the cells**

428 *Plate layout for the 8a T-Screen assay is presented in Figure 5.*

429 For every test item the half plate is prepared, as follows:

- 430 - Add 100 µL *working solutions* (*the double desired exposure concentrations*) of test item (TI) or
431 the appropriate controls prepared in PCM medium to the 100 µL already present (with the
432 cells) in the respective well(s) (to achieve a 1:1 dilution). Test all samples in triplicate with

- 433 exception of solvent control (SC) that is tested in 6-12 repetitions spread out each side of the
434 plate.
- 435 - Add 100 μ L 2x Ref(T3)C1 = the highest concentration in the range described in Section
436 3.2.2 Add 100 μ L PCM medium to three wells designed UC (PCM) (untreated control/cell
437 system control; Figure 5).
- 438 - Add 100 μ L of PCM medium to the 100 μ L already present (without the cells) in each well in
439 columns 1 and 12 (*Blank* and *Blank w/o AB*; Figure 5) that are used in the proliferation assay
440 for calculations.
- 441 - Add 100 μ L of PCM medium or PBS to each well in rows A and H (external wells) to ensure
442 proper humidity for the cells.
- 443 - Incubate plates for 96 ± 2 h at 37 ± 1 °C and $5 \pm 0.5\%$ (v/v) CO₂ in a humid atmosphere.

CONFIDENTIAL

444 **Figure 5.** Plate layout for the 8a T-Screen assay.

445 (A) the first plate is used to assess the range of Ref (T3) concentrations (upper part of the plate) and test item 1 (lower part of the plate) together with all
446 appropriate controls.

447

REF (T3)	Reference item T3 for AGONISM; Tested as full dose response curve once per series of plates. Each plate should contain the triplicate samples used in data analysis on the plate
PC (A)	Positive control; C max; Tested as triplicate sample once per series of plates
NC	Negative Control; C max; Tested as triplicate sample once per series of plates
UC (PCM)	Untreated Control; Tested as triplicate sample on each plate, PCM Medium Only,
SC	PCM medium + solvent; Tested as triplicate sample on each plate. It is recommended to use the same solvent for all reference and control items and to use that solvent for the test items, to ensure all samples are tested under the same conditions. In case another solvent is used for the test item, the impact on the test system and the results must be assessed.
SDS	Cytotox control; Tested as triplicate sample once per series of plates
*	Outer wells with PBS or PCM medium only; Outer wells may only be used only when the plate has the additional collar/space to fill it with liquid to ensure the cells proper humidity.
TI	Test item; Tested as full dose response curve
C	Concentrations of TI or REF (from C7/Cmin to C1/Cmax; C7 is the lowest concentration tested)
Blank	PCM medium without cells (AlamarBlue is added in the proliferation assay); used for calculations of % DR according to formulas given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.
Blank w/o AB	PCM medium without cells (AlamarBlue is not added in the proliferation assay); used for calculations of %DR according to formula given in SOP <i>Determination of cell proliferation in T-screen assay</i> in Section 2.

448

449 **Figure 5. continued**

450

451 **(A) the first plate:**

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
C	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
D	Blank	SC	Ref(T3) C7	Ref(T3) C6	Ref(T3) C5	Ref(T3) C4	Ref(T3) C3	Ref(T3) C2	Ref(T3) C1	UC (PCM)	SDS	Blank w/o AB
E	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
G	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
G	Blank	NC	TI 1 C7	TI 1 C6	TI 1 C5	TI 1 C4	TI 1 C3	TI 1 C2	TI 1 C1	PC(A)	SC	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

452

453

454 **(B) the next plates:**

455

	1	2	3	4	5	6	7	8	9	10	11	12
A	*	*	*	*	*	*	*	*	*	*	*	*
B	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
C	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
D	Blank	SC	TI 2 C7	TI 2 C6	TI 2 C5	TI 2 C4	TI 2 C3	TI 2 C2	TI 2 C1	UC (PCM)	SC	Blank w/o AB
E	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
G	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
G	Blank	SC	TI 3 C7	TI 3 C6	TI 3 C5	TI 3 C4	TI 3 C3	TI 3 C2	TI 3 C1	REF(T3) C1	SC	Blank w/o AB
H	*	*	*	*	*	*	*	*	*	*	*	*

456

457

458 **3.3.4 Measurement of cell proliferation**

459 Perform the cell proliferation assay according to Section 2: "AlamarBlue assay" in SOP
 460 "Determination of cell proliferation in T-screen assay". Before performing the assay, observe the cells
 461 under the microscope to record cytotoxic effect or precipitates (if any) as described in in SOP
 462 *Determination of cell proliferation in T-screen assay* in Section 2.1. Wells where precipitates were
 463 observed are excluded from calculation as recommended in SOP "Determination of cell
 464 *proliferation in T-screen assay*", Section 2.1.1 .

465 Observation of GH3 cell morphology after the exposure (focus on cytotoxic features, e.g. destroyed
 466 cells, destruction of the cell layers; intracytoplasmatic granules) is essential to confirm cytotoxic
 467 effect (the lower cellular activity of GH3 cells cultured in PCM medium without T3 as determined
 468 with the cell proliferation assay) and distinguish it from lack of agonistic effect (basal or low activity
 469 in PCM without T3).

470 **3.4 Data Analysis and Calculations**

471 Results of the cell proliferation assay, i.e. **% Dye reduction (%DR)**, calculated based on optical
 472 density (OD; absorbance) for each triplicate sample using formulas given in SOP *Determination of*
 473 *cell proliferation in T-screen assay* (Section 2.) are used to calculate the agonistic effect according to
 474 formulas given below (Section 3.4).

475 3.4.1 Based on received results (optical density/absorbance), calculate **% Dye reduction (%DR)** for
 476 each triplicate sample using formulas given in SOP *Determination of cell proliferation in T-*
 477 *screen assay* in Section 2.

478 3.4.2 Calculate the increase of cell proliferation that is expressed as **the relative proliferative**
 479 **effect (RPE)** for all samples. The response observed at 2 nM T3 [Ref(T3) C1] is considered as
 480 the maximum response and set as 100%. The response for the solvent control [SC] is set at
 481 0%. Exemplary results are presented in Figure 6A.

482
$$RPE = \frac{\%DR_x - \%DR_{SC}}{\%DR_{Ref(T3)C1} - \%DR_{SC}} \times 100$$

483 where:

- 484 SC - solvent control;
 485 Ref(T3) C1 - the max. response observed for T3 (the concentration no. 1; 2 nM);
 486 x - the effect of TI, PC, NC or REF T3 at the concentration analysed

487 the proliferation assay (SOP *Determination of cell proliferation in T-screen assay* in Section 2.)

488 3.4.3 Determine **the EC50 value of the reference item (Ref(T3))** and **the EC50 value of the test**
 489 **item (TI)** (for dose response) or the concentration giving highest RPE value above specific
 490 limit (e.g. > 10% increase or significantly different from SC). The EC50 value can be
 491 determined using the Hill curve model in a statistic programme e.g. GraphPad. The Hill curve
 492 model is a logistic regression model (variable slope, 4 parameters) that uses the following
 493 function:

494

$$y = RPE_{\min} + \frac{(RPE_{\max} - RPE_{\min})}{(1 + 10^{((\text{Log}EC_{50} - x) * \text{HillSlope}))}}$$

495

where:

496

x - Log of concentration

497

y - Relative induction (%)

498

RPE_{max} - Maximum relative induction (%)

499

RPE_{min} - Minimum relative induction (%)

500

LogEC₅₀ - Log of concentration at which 50% of maximum relative induction is observed

501

HillSlope - Slope factor of the Hill curve

502

3.4.4 Calculate the Z-factor for each plate tested according the following formula:

503

504

$$\mathbf{Z\text{-factor}}_{\text{plate no.}} = 1 - 3 \times \frac{(\text{SD}_{\text{plate no. [SC]}} + \text{SD}_{\text{plate no. [Ref(T3)C1]})}{\text{abs}(\%DR_{\text{plate no. [SC]}} - \%DR_{\text{plate no. [Ref(T3)C1]})}$$

505

where:

506

abs - absolute value;

507

SC - solvent control;

508

Ref (T3) C1 - the max. response observed for T3 (the concentration no. 1; 2 nM);

509

510


3.4.5 Relative potency of test items

511

The relative potency of test items is calculated by dividing the EC₅₀ of T3 by the EC₅₀ of the test item (or the concentration giving the highest RPE value that is above 10% increase as compared to SC).

513

Figure 6. Exemplary results obtained for reference item with agonistic potential



514

515

516

517 **3.5 Acceptance criteria**

518 **To be developed on basis of historical data with the reference and control items.**

519 3.5.1. Mean EC₅₀ value of T3 should be $-10 \pm 0.4 \log_{10}(\text{Molar})$ units (in the range from -10.4 to -9.6
520 $\log_{10}(\text{Molar})$ units).

521 3.5.2. %DR for TI SC and for Ref SC should not differ more than 15% from UC and REF SC

522 3.5.3. Z-factor > 0.5

523

524 **4. Reference(s)**

525 • Ghisari, M., & Bonefeld-Jorgensen, E. C. (2005). Impact of environmental chemicals on the
526 thyroid hormone function in pituitary rat GH3 cells. *Molecular and Cellular Endocrinology*, 244:
527 31–41. <http://doi.org/https://doi.org/10.1016/j.mce.2005.01.013>

528 • Gutleb, A. C., Meerts, I. A. T. M., Bergsma, J. H., Schriks, M., & Murk, A. J. (2005). T-Screen as a
529 tool to identify thyroid hormone receptor active compounds. *Environmental Toxicology and*
530 *Pharmacology*, 19: 231–238. <http://doi.org/https://doi.org/10.1016/j.etap.2004.06.003>

531 • Sirbasku, D.A., Pakala, R., Sato, H., Eby, J.E. (1991). Thyroid hormone dependent pituitary tumor
532 cell growth in serum-free chemically defined culture. A new regulatory role for apotransferrin.
533 *Biochem.* 30: 7466–7477.

534 • Schriks M, Vrabie CM, Gutleb AC, Faassen EJ, Rietjens IM, Murk AJ. (2006) T-screen to quantify
535 functional potentiating, antagonistic and thyroid hormone-like activities of poly halogenated
536 aromatic hydrocarbons (PHAHs). *Toxicol In Vitro.* 20: 490-498.

537 • Hohenwarter O, Waltenberger A, Katinger H. (1990) An in vitro test system for thyroid hormone
538 action. *Anal. Biochem.* 234: 56-59.

539 • Koyanagi M, Kawakabe S, Arimura Y. (2016). A comparative study of colorimetric cell proliferation
540 assays in immune cells. *Cytotechnology*, 68, 1489-1498. doi:10.1007/s10616-015-9909-2.

541

542 **5. A change log**

543

544 Standard Operating Procedure:

545 T-screen assay measuring cell proliferation of GH3 cells using alamar blue/resazurin

546

ID	Version	Status	Date	List of changes in version 02 with respect to version 01:
1.	v. 01	Approved	03.03.2021	not applicable
2.	v. 02	Draft	09.12.2021	- Antagonist part of the method was removed throughout the SOP - %AR calculation was removed as the results were the same as for the %DR calculation
3.	v. 02	Draft	27.12.2021	- Adding a section: A change log
4.	v. 02	Final	25.01.2022	- updated incubation time from 48h ± 1h to 48 ± 2h, and updated exposure time from 96 ± 1h to 96 ± 2h as it is not critical and allows for more flexibility'

547

CONFIDENTIAL

GETTING IN TOUCH WITH THE EU

In person

All over the European Union there are hundreds of Europe Direct centres. You can find the address of the centre nearest you online (european-union.europa.eu/contact-eu/meet-us_en).

On the phone or in writing

Europe Direct is a service that answers your questions about the European Union. You can contact this service:

- by freephone: 00 800 6 7 8 9 10 11 (certain operators may charge for these calls),
- at the following standard number: +32 22999696,
- via the following form: european-union.europa.eu/contact-eu/write-us_en.

FINDING INFORMATION ABOUT THE EU

Online

Information about the European Union in all the official languages of the EU is available on the Europa website (european-union.europa.eu).

EU publications

You can view or order EU publications at op.europa.eu/en/publications. Multiple copies of free publications can be obtained by contacting Europe Direct or your local documentation centre (european-union.europa.eu/contact-eu/meet-us_en).

EU law and related documents

For access to legal information from the EU, including all EU law since 1951 in all the official language versions, go to EUR-Lex (eur-lex.europa.eu).

Open data from the EU

The portal data.europa.eu provides access to open datasets from the EU institutions, bodies and agencies. These can be downloaded and reused for free, for both commercial and non-commercial purposes. The portal also provides access to a wealth of datasets from European countries.

The European Commission's science and knowledge service

Joint Research Centre

JRC Mission

As the science and knowledge service of the European Commission, the Joint Research Centre's mission is to support EU policies with independent evidence throughout the whole policy cycle.



EU Science Hub
joint-research-centre.ec.europa.eu



@EU_ScienceHub



EU Science Hub - Joint Research Centre



EU Science, Research and Innovation



EU Science Hub



EU Science